oceanus

35

EN	BALLOON CATHETER FOR PERIPHERAL PREDILATION FOR 0.035" GUIDE WIRE. INSTRUCTIONS FOR USE	3
ES	CATÉTER BALÓN DE PREDILATACIÓN PERIFÉRICO PARA GUÍA DE ALAMBRE DE 0,035". INSTRUCCIONES DE USO	7
PT	CATETER BALÃO DE PRÉ-DILATAÇÃO PERIFÉRICO PARA FIO-GUIA DE 0,035". INSTRUÇÕES DE UTILIZAÇÃO	11
IT	CATETERE A PALLONCINO PERIFERICO DI PRE-DILATAZIONE PER FILO GUIDA DI 0,035".ISTRUZIONI D'USO	15
FR	CATHÉTER BALLONNET DE PRÉ-DILATATION PÉRIPHÉRIQUE POUR GUIDE EN FIL MÉTALLIQUE DE 0,035" MODE D'EMPLOI	19
DE	BALLONKATHETER FÜR PERIPHERE VORDILATATION FÜR 0,035" FÜHRUNGSdraht. GERBRAUCHSHINWEISE	23
RU	БАЛЛОНЫЙ КАТЕТЕР ДЛЯ ПРЕДИЛАТАЦИИ ПЕРИФЕРИЧЕСКИХ СОСУДОВ ДЛЯ ПРОВОДНИКА 0,035". ИНСТРУКЦИИ ПО ИСПОЛЬЗОВАНИЮ	27
CS	BALÓNKOVÝ KATÉTR PRO PERIFERNÍ PREDILATACI S VODICÍM DRÁTEM 0,035". NÁVOD K POUŽITÍ	31
EL	ΚΑΘΕΤΗΡΑΣ ΜΕ ΜΠΑΛΟΝΙ ΓΙΑ ΠΡΟΔΙΑΣΤΟΛΗ ΠΕΡΙΦΕΡΙΚΩΝ ΑΡΤΗΡΙΩΝ ΓΙΑ ΣΥΡΜΑΤΙΝΟ ΟΔΗΓΟ 0,035 τι. ΟΔΗΓΙΕΣ ΧΡΗΣΗΣ	35
PL	CEWNIK BALONOWY DO PREDYLATACJI NACZYŃ OBWODOWYCH DO PROWADNIKÓW 0,035 cala. INSTRUKCJA UŻYTKOWANIA	39

INSTRUCTIONS FOR USE

BALLOON CATHETER FOR PERIPHERAL PREDILATION FOR 0.035" GUIDE WIRE.

1. Description of the Product

The balloon catheter is a dual lumen catheter from the connector to the tip, also known as OTW, designed for percutaneous transluminal angioplasty of peripheral arteries.

The catheter has two lumens body from the connector to the tip.

On the distal part of the catheter, before the tip, is located the balloon or inflatable segment to different pressures. The balloon dilates the artery when is inflated by means of infusion of contrast liquid inside it.

The connector is Y-shaped, with two ports:

- The side port allows contrast liquid passage to dilate the balloon.
- The straight port intended for guide wire insertion. The guide wire is inserted into the connector port and comes out of the tip of the catheter to reach the lesion.

There is likewise a lumen through which the guide wire runs and a lumen through which contrast liquid enters the balloon, directly connected to each of the ports in the luer or connector.

The maximum diameter of the guide wire must not exceed 0.89mm = 0.035 inches.

The useful length of the catheter is 80cm or 140cm depending on the item reference.

The balloon is supplied in the following lengths and diameters:

		BALLOON DIAMETERS (mm)								
		3	4	5	6	7	8	9	10	12
BALLOON LENGTH (mm)	20	X	X	X	X	X	X	X	X	X
	40	X	X	X	X	X	X	X	X	X
	60	X	X	X	X	X	X	X	X	X
	80	X	X	X	X	X	X	X		
	120	X	X	X	X	X				
	150	X	X	X	X					
	200	X	X	X	X					

The balloon is intended to reach different diameters at different pressures (according to the compliance curve shown on the primary packaging).

To inflate the balloon, the connector inflation port must be connected to an inflation device. The balloon is inflated to predictable diameters with the pressure shown by the manometer. Two radiopaque markers are located at each end of the balloon in order to mark its length and help the user to see the catheter while inside the patient.

At the distal end is the tip, which is rounded andatraumatic in shape in order to avoid damaging the arteries while it is being advanced.

Contents

- One OTW catheter covered by protection sheath for the balloon, inserted into a dispenser in a sterile bag.
- One card with the compliance curve showing the nominal inflation pressure and the recommended maximum pressure.
- One leaflet with instructions for use.

2. Indications

Indicated for dilation of stenosis located in the renal, iliac, femoral, popliteal and infrapopliteal arteries, and for treating obstructive lesions of native or synthetic arteriovenous dialysis fistulae.

3. Contraindications

Inability to cross lesion with a guide wire.

INSTRUCTIONS FOR USE**BALLOON CATHETER FOR PERIPHERAL PREDILATION FOR 0.035" GUIDE WIRE.****4. Warnings**

- The device should only be used by experienced doctors in the percutaneous transluminal angioplasty (PTA) technique.
- The percutaneous procedure must be carried out in hospitals equipped with emergency facilities to conduct open heart surgery or hospitals with easy access to other hospitals where this type of operation can be performed.
- The product is for single use. It must not be re-sterilised or re-used once the procedure has been completed. Re-using the product in another patient may lead to cross contamination, infection or transmission of infectious diseases from one patient to another. Re-use of the product may cause alterations to it and affect its effectiveness.
- The product is supplied sterile. Check the expiry date and do not use products which are past this date.
- Inspect the packaging before opening it. If there are any defects or the packaging is damaged, do not use the product.
- Use aseptic techniques when the product is taken out of the primary packaging.
- Select the balloon diameter according to the diameter of the artery to be treated. Do not insert a balloon with a diameter larger than the artery to be treated.
- Do not dry with gauze.
- Do not expose the product to organic solvents or to contrast media not suitable for intravascular use.
- Do not use air or gaseous media to inflate the balloon.
- Advance the product on the guide wire using fluoroscopy. Do not allow the product to advance without the guide wire inside it.
- Do not manipulate, advance or retract either the catheter or the wire guide while the balloon is inflated.

5. Precautions

- Administer suitable medical therapy to the patient: anticoagulants, vasodilators, etc., in accordance with the procedure for inserting intravascular catheters.
- Before use, check the size and suitability of the catheter for its intended use.
- Proceed with extreme care so as not to damage the catheter while advancing it.
- The balloon must be inflated with a mixture of saline solution and contrast liquid (preferably 50/50, though the proportion of contrast liquid can be reduced in large balloons).
- Do not exceed the rated burst pressure (RBP) as the balloon could burst.
- If any resistance is noted during advance, stop and determine its cause before continuing.
- If any resistance to withdrawing the catheter is noted, it is recommended that the balloon catheter, the guide wire and introducer be removed as a single unit.
- Store in a cool, dry place away from direct sunlight.
- After use, this product may represent a biohazard. Handle and dispose it in accordance with accepted medical practices and pertinent local, state or federal laws and regulations.

6. Possible Adverse Effect / Complications

The possible adverse effects and/or complications which might arise before, during or after the procedure include the following:

- Death
- Ictus / embolism / thrombosis
- Haemodynamic deterioration
- Acute elastic recoil
- Arrhythmia

INSTRUCTIONS FOR USE**BALLOON CATHETER FOR PERIPHERAL PREDILATION FOR 0.035" GUIDE WIRE.**

- Aneurysm or pseudo-aneurysm
- Infections
- Total occlusion of the artery
- Reocclusion of the treated area: restenosis
- Spasm
- Perforation or dissection of the treated area
- Arteriovenous fistula
- Local bleeding with haematoma at the access site
- Hypo/hypertension
- Allergic reactions to materials.

7. How to Use**7.1. Equipment Required**

- Heparinized normal saline solution.
- Contrast medium (use contrast media suitable for intravascular use).
- 0.035" guide wire; do not use any other size of guide wire.
- Introducer with a haemostatic valve of the size specified on the label. Do not use smaller sizes as these might damage the catheter. If an introducer more than 25cm long or of the braided type is used, it may be necessary to increase its size to avoid rubbing on the catheter.
- Three-way stopcock.
- Inflation device (manual pump with built-in manometer).
- Several standard 10-20cc syringes with saline solution for washing the system.

7.2. Preparing the Catheter

- Remove the catheter from the protective dispenser. Check that it is the right size. Remove the protective sheath from the balloon.
- Wash the product with a sterile isotonic saline solution or similar on the outside and inside the channel through which the guide wire passes, before inserting it into the patient.
- Connect a 10ml syringe containing sterile saline solution to the port for the guide wire (straight part of the connector) and irrigate the lumen until the liquid comes out of the tip.
- Purge the whole system of air:
 - Attach a three-way stopcock to the contrast liquid connector port (the side part of the connector).
 - Close the air passage through the balloon.
 - Attach a 10-20ml syringe with a third part of saline solution to the three-way stopcock.
 - Open the passage in the three-way stopcock between the syringe and the catheter.
 - With the syringe in a vertical position, withdraw the plunger upwards allowing the air bubbles to exit to the liquid.
 - When bubbles stop coming in, close the three-way stopcock at the catheter end and remove the syringe.

DANGER: If bubbles do not stop coming from the catheter into the syringe after 3 minutes of negative pressure, this may be a clear indication that the balloon catheter has leaks, is broken or the connections between the syringe and the three-way key are not properly sealed. If bubbles continue to be seen after checking the connections, do not use the device. Return it to the manufacturer or distributor for inspection.

7.3. Technique for Insertion / Treatment

- Position the 0.035" guide wire across the lesion using PTA techniques, using fluoroscopy to determine its position at all times.

**INSTRUCTIONS FOR USE
BALLOON CATHETER FOR PERIPHERAL PREDILATION FOR 0.035" GUIDE WIRE.**

- Insert the proximal end of the guide into the distal tip of the balloon catheter.
- Carefully advance the balloon catheter on the guide wire, through the introducer and with the aid of fluoroscopy, until it reaches the area to be treated. In the event that the balloon catheter selected cannot be placed inside the lesion, take a smaller size to predilate the lesion.
- Connect the inflation device (1/3 of its capacity of contract liquid mixture: saline solution) to the three-way stopcock and open the passage between the balloon and the inflation device.
- Apply pressure to the inflation device so that the balloon inflates. Do not exceed the rated burst pressure (RBP) on the label and the compliance curve.
- Apply pressure for 30-60 seconds depending on the lesion.
- Draw the plunger on the inflation device back to deflate the balloon. Keep up the negative pressure for 20 to 60 seconds depending on the size of the balloon. Make sure the balloon is fully deflated (by means of fluoroscopy) before moving the catheter.
- With negative pressure in the inflation device and with the guide wire in position, carefully retract the catheter until it can be removed through the introducer.
- Remove the guide wire.

8. Warranty

The device and all its components have been designed, manufactured, tested and packed with a maximum level of precaution. LVD Biotech warrants the device until its expiration date, provided the packaging is not broken, tampered with or damaged.

**STANDARDISED SYMBOLS FOR USE IN LABELLING MEDICAL DEVICES
(EN 980:2008 STANDARD)**

Catalogue reference

Batch no.

Expiry date

Keep away from sunlight

Manufacturer

Sterilised with ethylene oxide

Limitation in storage temperature

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© iVascular - cod. MP45001 ed2 Jun 2015