luminor 35

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**INSTRUCTIONS FOR USE. PACLITAXEL ELUTING
PERIPHERAL BALLOON DILATATION CATHETER, FOR 0.035" GUIDEWIRE****1. Description of the Product**

The Paclitaxel eluting dilatation balloon is a double lumen catheter from the connector to the tip, also called OTW, designed for percutaneous transluminal angioplasties of peripheral arteries.

The catheter has two lumens body from the connector to the tip.

On the catheter's distal part, before the tip, there is the balloon with the drug or segment inflatable at different pressures that dilates the artery upon inflating by means of infusion of contrast fluid inside it. The connector is Y-shaped and it has two entry ports.

The lateral port, which connects with the smallest lumen of the catheter's body, allows entry of contrast fluid to dilate the balloon.

The straight port is the port for guidewire passage. The guidewire is introduced in the connector's entry port and going through the largest lumen, it goes out through the catheter's tip to take it to the lesion.

There is likewise a lumen through which the guide wire runs and a lumen through which contrast liquid enters the balloon, directly connected to each of the ports in the luer or connector.

The maximum diameter of the guide wire must not exceed 0.89mm = 0.035 inches.

The balloon is coated with the drug Paclitaxel, included in a physiologically innocuous matrix. The drug's dose is 3 micrograms / mm² of balloon surface and it is intended to avoid cellular proliferation, consequently decreasing the re-intervention rate.

The drug is released from the balloon by means of a rapid inflation in a way so that a high dose of the drug is released in a very short period of time. For releasing the drug onto the artery wall, the inflation process must last from 30 seconds to 1 minute, but dilation of the lesion must be optimized by using longer inflation times at the operator's discretion.

The useful length of the catheter is 80cm or 140cm depending on the item reference.

The balloon is supplied in the following lengths and diameters:

		BALLOON DIAMETERS (mm)		
		5	6	7
BALLOON LENGTH (mm)	20	X	X	X
	40	X	X	X
	60	X	X	X
	80	X	X	X
	120	X	X	X
	150	X	X	

The balloon is intended to reach different diameters at different pressures (according to the compliance curve shown on the primary packaging).

To inflate the balloon, it is necessary to connect the connector's inflation port to a syringe that includes a manometer. With the pressure indicated by the manometer, the balloon is inflated to predictable diameters.

At both ends of the balloon, two radiopaque markers are included so as to delimit its length, and help the user know the place where it is when it is inside the patient.

At the distal end is the tip, which is rounded and atraumatic in shape in order to avoid damaging the arteries while it is being advanced.

Contents

- A paclitaxel releasing dilatation balloon catheter OTW, covered by a protector over the balloon and introduced in a dispenser and packed in a sterile bag.
- One card with the compliance curve showing the nominal inflation pressure and the recommended maximum pressure.
- One leaflet with instructions for use.

**INSTRUCTIONS FOR USE. PACLITAXEL ELUTING
PERIPHERAL BALLOON DILATATION CATHETER, FOR 0.035" GUIDEWIRE**

EN

2. Indications

The device is indicated for dilation of stenosis located in the iliac, femoral, iliofemoral, popliteal, infrapopliteal and renal arteries, as well as for the treatment of obstructive lesions of arteriovenous fistulas, whether original or artificial, with a reference diameter balloon from 5 to 7 mm and lengths from 20 to 150 mm. It is also indicated for stent postdilation in the peripheral vascular system.

3. Contraindications

Inability to cross lesion with a guide wire.

It must not be used in pregnant women or during lactation or in patients with known hypersensitivity to paclitaxel.

4. Warnings

- The device should only be used by experienced doctors in the percutaneous transluminal angioplasty (PTA) technique.
- The percutaneous procedure must be carried out in hospitals equipped with emergency facilities to conduct open heart surgery or hospitals with easy access to other hospitals where this type of operation can be performed.
- The product is for single use. It must not be re-sterilised or re-used once the procedure has been completed. Re-using the product in another patient may lead to cross contamination, infection or transmission of infectious diseases from one patient to another. Re-use of the product may cause alterations to it and affect its effectiveness.
- The product is supplied sterile. Check the expiry date and do not use products which are past this date.
- Inspect the packaging before opening it. If there are any defects or the packaging is damaged, do not use the product.
- Use aseptic techniques when the product is taken out of the primary packaging.
- Select the balloon diameter according to the diameter of the artery to be treated. Do not insert a balloon with a diameter larger than the artery to be treated.
- Do not touch the balloon directly or dry it with gauze, since doing so may cause the drug coating to come off.
- Do not expose the product to organic solvents such as alcohol, or to contrast media not indicated for intravascular use.
- Do not use air or gaseous media to inflate the balloon.
- Advance the product on the guide wire using fluoroscopy. Do not allow the product to advance without the guide wire inside it.
- Do not manipulate, advance or retract either the catheter or the wire guide while the balloon is inflated.
- In case a post-dilatation is necessary after use of product, it must be made with a conventional dilatation balloon, one without the drug, since the same vessel segment must not be dilated with more than one balloon with drug.
- Drug eluting stents must not be implanted in the same segment that has been previously treated with the balloon with drug, since it is possible for overdose or interaction among the active principles to occur.
- Published in vitro and in vivo studies prove the mutagenic effect of Paclitaxel. The carcinogenesis risk potential has not been studied. The published in vivo studies prove the teratogenic effect, therefore it is not recommended for pregnant women to use this product.

5. Precautions

- Administer suitable medical therapy to the patient: anticoagulants, vasodilators, etc., in accordance with the procedure for inserting intravascular catheters.
- Before use, check the size and suitability of the catheter for its intended use.
- Use the product with precautions in complex lesions, since the abrasive nature of said lesions may alter the drug coating. In these cases it is recommended to predilate the lesion previously with a con-

**INSTRUCTIONS FOR USE. PACLITAXEL ELUTING
PERIPHERAL BALLOON DILATATION CATHETER, FOR 0.035" GUIDEWIRE**

ventional, drugless balloon.

- Extreme care must be exercised when the device is taken off from the pack and passed through the hemostatic valve to make sure that the drug coating is not damaged.
- Proceed with extreme care so as not to damage the catheter while advancing it.
- The balloon must be inflated with a mixture of saline solution and contrast liquid (preferably 50/50, though the proportion of contrast liquid can be reduced in large balloons).
- Do not exceed the rated burst pressure (RBP) as the balloon could burst.
- If any resistance is noted during advance, stop and determine its cause before continuing.
- If any resistance to withdrawing the catheter is noted, it is recommended that the balloon catheter, the guide wire and introducer be removed as a single unit.
- Store in a cool, dry place away from direct sunlight.
- After use, this product may represent a biohazard. Handle and dispose it in accordance with accepted medical practices and pertinent local, state or federal laws and regulations.

6. Possible Adverse Effect / Complications

The possible adverse effects and/or complications which might arise before, during or after the procedure include the following:

- Death
- Ictus / embolism / thrombosis
- Haemodynamic deterioration
- Acute elastic recoil
- Arrhythmia
- Aneurysm or pseudoaneurysm
- Infections
- Total occlusion of the artery
- Reocclusion of the area treated: restenosis
- Spasm
- Perforation or dissection of the area treated
- Arteriovenous fistula
- Local bleeding with haematoma at the access site
- Hypo/hypertension
- Allergic or immune reaction to the drug.
- Histological changes in the arterial wall.

7. Pharmacological interaction

The potential interactions of paclitaxel with the drugs administered concomitantly have not been investigated. Drug interactions at the systemic levels of paclitaxel are detailed in the product data sheet but, considering that the amount of paclitaxel on the balloon surface corresponds approximately to a few hundredths of the amount normally used in the antineoplastic treatment, it is considered very unlikely for interactions with other drugs to occur.

Nevertheless, caution must be exercised when CYP3A4 or CYP2C8 substrates (such as terfenadine, cyclosporin, lovastatin, midazolam and ondansetron) or drugs with high PPB (Plasma Protein Binding) are administered concomitantly.

8. How to Use**8.1 Equipment Required**

- Heparinized normal saline solution
- Contrast medium (use contrast media suitable for intravascular use)
- 0.035" guide wire: do not use any other size of guide wire
- Introducer with a haemostatic valve of the size specified on the label. Do not use smaller sizes as

**INSTRUCTIONS FOR USE. PACLITAXEL ELUTING
PERIPHERAL BALLOON DILATATION CATHETER, FOR 0.035" GUIDEWIRE**

these might damage the catheter. If an introducer more than 25cm long or of the braided type is used, it may be necessary to increase its size to avoid rubbing on the catheter.

- Three-way stopcock
- Inflation device (manual pump with built-in manometer)
- Several standard 10-20cc syringes with saline solution for washing the system

8.2 Preparing the Catheter

- In case of occlusion or very calcified lesions, it is recommended to predilate the lesion with a conventional balloon, without drug.
- Remove the catheter from the protection dispenser. Check that it is the right size. Remove the balloon protection sheath smoothly.
- Connect a 10ml syringe containing sterile saline solution to the port for the guide wire (straight part of the connector) and irrigate the lumen until the liquid comes out of the tip.
- Purge the whole system of air.
 - Attach a three-way stopcock to the contrast liquid connector port (the side part of the connector).
 - Close the air passage through the balloon.
 - Attach a 10-20ml syringe with a third part of saline solution to the three-way stopcock.
 - Open the passage in the three-way stopcock between the syringe and the catheter.
 - With the syringe in a vertical position, withdraw the plunger upwards allowing the air bubbles to exit to the liquid.
 - When bubbles stop coming in, close the three-way stopcock at the catheter end and remove the syringe.

DANGER: If bubbles do not stop coming from the catheter into the syringe after 3 minutes of negative pressure, this may be a clear indication that the balloon catheter has leaks, is broken or the connections between the syringe and the three-way key are not properly sealed. If bubbles continue to be seen after checking the connections, do not use the device. Return it to the manufacturer or distributor for inspection.

8.3 Technique for Insertion / Treatment

- Position the 0.035 guide wire across the lesion using PTA techniques, using fluoroscopy to determine its position at all times.
- Insert the proximal end of the guide into the distal tip of the balloon catheter.
- Advance carefully the balloon catheter along the guidewire, through the guide catheter and with the help of fluoroscopy, until the area to be treated is reached. Extreme care must be exercised when the balloon catheter is passed through the hemostatic valve to make sure that the drug coating is not damaged. In case the balloon catheter chosen cannot be placed inside the lesion, choose a smaller size of a drugless conventional balloon in order to predilate the lesion.
- Connect the inflation device (1/3 of its capacity of contract liquid mixture: saline solution) to the three-way stopcock and open the passage between the balloon and the inflation device.
- Apply pressure to the inflation device so that the balloon inflates. Do not exceed the rated burst pressure (RBP) on the label and the compliance curve.
- Keep the pressure for 30-60 seconds for an optimal release of the drug. Most of the drug is released in the first 30 seconds of inflation but, in order to optimize the dilation of the lesion, it is possible to use longer inflation times greater than 1 minute at the operator's discretion.
- Draw the plunger on the inflation device back to deflate the balloon. Keep up the negative pressure for 20 to 60 seconds depending on the size of the balloon. Make sure the balloon is fully deflated (by means of fluoroscopy) before moving the catheter.
- With negative pressure in the inflation device and with the guide wire in position, carefully retract

**INSTRUCTIONS FOR USE. PACLITAXEL ELUTING
PERIPHERAL BALLOON DILATATION CATHETER, FOR 0.035" GUIDEWIRE**

the catheter until it can be removed through the introducer.

- Remove the guide wire.
- Double antiplatelet therapy must be administered for a minimum of four weeks after the intervention, and the treatment may be prolonged at the operator's discretion.

9. Warranty

The device and all its components have been designed, manufactured, tested and packed with a maximum level of precaution. LVD Biotech guarantees the product up to its expiry date providing its packaging is not broken, tampered with or damaged.

**STANDARDIZED SYMBOLS FOR USE IN LABELLING MEDICAL DEVICES
(EN 980:2008 STANDARD)**

Catalogue reference

Batch no.

Expiry date

Keep away from sunlight

Manufacturer

Sterilized with ethylene oxide

Limitation in storage temperature

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