
# iVascular® angiolute BTK

EN	INSTRUCTIONS FOR USE SIROLIMUS-ELUTING PERIPHERAL STENT SYSTEM	2
ES	INSTRUCCIONES DE USO. SISTEMA DE STENT PERIFÉRICO CON LIBERACIÓN DE SIROLIMUS	4
PT	INSTRUÇÕES DE USO. SISTEMA DE STENT PERIFÉRICO COM LIBERAÇÃO DE SIROLIMUS	6
IT	ISTRUZIONI PER L'USO. SISTEMA DI STENT PERIFERICO A RILASCIO DI SIROLIMUS	8
FR	INSTRUCTIONS D'UTILISATION. SYSTÈME D'ENDOPROTHÈSE PÉRIPHÉRIQUE AVEC ELUTION DE SIROLIMUS	10
DE	GEBRAUCHSHINWEISE. SIROLIMUS-FREISETZENDES PERIPHERES STENTSYSTEM	12
RU	ИНСТРУКЦИИ ПО ПРИМЕНЕНИЮ. ПЕРИФЕРИЧЕСКАЯ СТЕНТОВАЯ СИСТЕМА С СИРОЛИМУСОВЫМ ПОКРЫТИЕМ	14
CS	NÁVOD K POUŽÍTÍ. SYSTÉM PERIFERNÍHO STENTU UVOLŇUJÍCÍHO SIROLIMUS	17
EL	ΟΔΗΓΙΕΣ ΧΡΗΣΗΣ. ΣΥΣΤΗΜΑ ΠΕΡΙΦΕΡΙΚΟΥ ΣΤΕΝΤ ΠΕΡΙΦΕΡΕΙΑΚΟΣ ΑΓΓΕΙΩΝ ΜΕ ΑΠΕΛΕΥΘΕΡΩΣΗ SIROLIMUS	19
PL	NSTRUKCJA OBSŁUGI. SYSTEM STENTU OBWODOWEGO UWALNIAJĄCEGO SIROLIMUS	21
DU	GEBRUIKSAANWIJZING SIROLIMUS-ELUEREND PERIFEER STENTSYSTEEM	23
SE	BRUKSANVISNING TILL SIROLIMUS-ELUERANDE PERIFERT STENTSYSTEM	25
TU	KULLANMA TALIMATI SİROLİMUS SALINIMLI PERİFERİK STENT SİSTEMİ	27

## ENGLISH

### 1. Description of the Product

The LVD Biotech SL angiolite BTK Sirolimus-eluting peripheral stent is a stent made of a Chromium cobalt alloy named L605, to which a coating made of a mixture of sirolimus and polymers has been applied. The stent is supplied pre-mounted on a delivery system which will expand it, by means of the distal balloon, in the artery to be treated.

The stent is intended to fit with different artery diameters by adjusting its open-cell design with alternated connecting bridges. The stent is manufactured means of laser cuttings in metallic tubes. Afterwards, it is submitted to several treatments in order to give a smooth and sleek finishing to the surface.

The content of a stent drug is approximately 1.4 µg/mm<sup>2</sup> with the stent being larger (4.50x39 mm), which has the maximum drug amount, 339 µg.

The stent delivery system consists on a rapid-exchange balloon catheter. The catheter comprises an inflatable balloon at the distal part. The balloon has different diameters and lengths to fit with the different stent configurations. When the balloon is inflated, it dilates the stent and deploys it against the artery. Afterwards, the balloon is deflated and removed, and the stent remains permanently implanted.

There are markers on the catheter shaft to help the user to calculate the position of the catheter as it is advanced through the guiding catheter without fluoroscopy, so that, when the last mark disappears, the catheter is already near to the tip of the guiding catheter and about to enter the artery. The marker located nearest to the catheter adapter is for femoral guiding catheters and the one furthest away is for brachial guiding catheters.

The distal part of the catheter is covered in a durable hydrophilic coating to lubricate the catheter so that it can navigate through arteries.

The maximum guidewire passage diameter must not exceed 0.36 mm = 0.014".

Useful length of the catheter is 142 cm, while the total length is 150 cm.

The stent is supplied in the following stent lengths and diameters: (pg. 31).

To expand the stent, the proximal luer-lock hub is needed to be connected to an inflation device with a manometer. The balloon is inflated with the pressure shown by the manometer to expand the stent to predictable diameters. One radiopaque marker is located at each end of the balloon in order to mark its length and help the user to know its position while it is inside the patient.

### Contents

- One sirolimus-eluting peripheral stent system comprising the coated stent pre-mounted on its delivery system (balloon catheter). The stent is protected by a sheath and a stylet inserted at the guide insertion lumen, and it is inserted into a dispenser and in a sterile bag. In turn, the bag with sterile contents is inserted into a second foil bag.
- One card with the compliance curve showing the working pressure range
- One implant card
- One leaflet with instructions for use

### 2. Indications

The device is indicated for treating chronic and acute arterial lesions in the lower limbs below the knee (BTK), as well as popliteal and infrapopliteal lesions with reference vessel diameters of between ≥2 mm and ≤4.50 mm, with the aim of increasing the internal diameter of an artery and consequently improving blood flow in the following cases:

- Severe claudication
- Critical limb ischaemia

### 3. Contraindications

- Antiplatelet and/or anticoagulation therapy is contraindicated.
- Patients with clear diagnosed heavy metal allergy.
- Patients with lesions that prevent complete angioplasty balloon inflation or proper stenting.
- Patients with hypersensitivity or an allergy to Sirolimus or its by-products and to acrylate polymers.

### 4. Warnings

- Stent implantation must be carried out in hospitals equipped with emergency facilities to conduct open heart surgery or, in its absence, in

hospitals with easy access to other hospitals where this type of operation can be performed.

- The product should only be used by experienced physicians in hemodynamic and familiar with the percutaneous peripheral intervention and stent implantation.
- The device is designed for single use in a single patient. DO NOT RE-STERILIZE OR RE-USE IT, REPROCESSING IS FORBIDDEN. Re-using the product in another patient may lead to cross contamination, infection or transmission of infectious diseases from one patient to another. Re-use of the product may cause alterations to it and affect its effectiveness.
- The product is supplied sterile. Check its expiry date and do not use products which are past this date.
- Administer suitable medical therapy to the patient: anticoagulants, vasodilators, etc., according to the protocol of peripheral intravascular catheter insertion.
- Use aseptic techniques when the product is taken out of its primary packaging.
- Do not dry with gauze.
- Do not expose the delivery device to organic solvents.
- Do not use oily or high viscous contrast media unsuitable for intravascular use.
- Do not use air or gaseous media to inflate the balloon, it must be inflated with a mixture of saline solution and contrast liquid (preferably 1:1).
- Select the suitable product size, in terms of diameter and length, according to the size of the lesion fluoroscopically observed.
- Advance the product on the guidewire by means of fluoroscopy imaging. Do not allow product advance without the guidewire inside it.
- Do not overlap stents with different compositions.
- The possible interaction of the product with other drug-eluting stents has not been evaluated, and therefore, they should be avoided whenever possible.
- Do not expose the product to hospital waste material.
- Studies in pregnant or breastfeeding women or in men who wish to have children have not been conducted with this product. Fertility studies in animals demonstrated reproductive toxicity. The use of the product in women who are trying to become pregnant or who are pregnant or breastfeeding is not recommended.

### 5. Precautions

#### 5.1 General precautions

- It is recommended to be stored in a cool, dry place (15° -30° C), away from direct sunlight.
- Inspect the packaging. If it is damaged, do not use the product.
- After its use, this product may be a biological hazard. Handle and dispose it of according to the accepted medical practices, and the pertinent local, state or federal laws and standards.
- The foil bag is not sterile. The sterile protective barrier is the plastic bag.
- Safety and efficacy of the product in patients whose lesion has been previously treated with brachytherapy have not yet been established.
- The use of drug-eluting stents in patients and lesions not included in the instructions for use may mean a higher risk of adverse reactions, such as stent thrombosis, myocardial infarction or death.
- The use of drug-eluting stents in patients and lesions not listed in the instructions for use may lead to an increased risk of adverse reactions.

#### 5.2 Precautions during product preparation

- Inspect the product before using it to ensure it is in perfect conditions, and discard any device not showing proper confidence.
- Do not prepare or pre-inflate the delivery system before starting the procedure or deploying the stent at the lesion.
- Do not handle the stent with the fingers, due to this action may dislodge the stent from the balloon and damage the coating.
- Do not attempt to separate the stent from its delivery system; if any stent movement or unsuitable stent place or crimping is noticed, do not use it.

- Do not attempt to straighten a kinked hypotube (proximal part of the delivery system), it may get broken.

### 5.3 Precautions during product advance-positioning-stenting-withdrawal

- In case of multiple lesions treatment, begin with the most distal lesions, followed by the proximal ones.
- Always advance the device with fluoroscopic imaging. If any resistance is noted, do not continue until determining its cause. If the stent is not able to cross or reach the lesion, remove the whole system (including the guiding catheter and the guidewire) as a single unit.
- If any resistance is noticed during advance, stop the advance and determine its cause before continuing.
- Once stent expansion is started, do not intent to remove it or change its position.
- Do not manipulate, advance or retract, either the catheter or the guidewire, while the balloon is inflated.
- Do not exceed the maximum recommended pressure indicated on the label and the attached compliance curve.
- If any resistance is noted during lesion access, delivery system withdrawal or withdrawal of a stent which could not be implanted, the whole system should be removed as a single unit:
- Do not withdraw the delivery system into the guiding catheter or long introducer:
  - Place the proximal balloon marker distally to the tip of the guiding catheter or long introducer
  - Advance the guidewire towards into the anatomy as far as safely possible
  - Tighten the haemostatic valve so it properly grabs the guiding catheter or long introducer and the balloon
  - Remove the guiding catheter or long introducer and the balloon catheter together, as a single unit
  - Finally, remove the guidewire, or in case an implantation is wished to begin again, maintain the guidewire and placed the rest of the devices.
- Great care must be exercised immediately after implantation and if it was needed to re-cross it, not to push the stent with guidewires or balloons not to destructure it.

### 6. Precautions: MRI[1] Safety:

Chromium-cobalt drug-eluting stents have shown in pre-clinical trials, individually and in super-posed configurations, that they are safe for magnetic resonance under certain conditions:

- Static magnetic field with intensity equal to or less than 3 Tesla strength
- Maximum spatial gradient of 720 gauss/cm
- Whole-body specific absorption rate (SAR) with a maximum average of 3.0 W/Kg over 15 minutes of magnetic resonance exploration.

In non-clinical tests with static magnetic resonance of 3.0 Tesla strength and a whole-body maximum specific absorption rate of 2.0 W/Kg over 15 minutes of magnetic resonance exploration, the increase in temperature for a single stent was a little more than 1°C, while for an overlapping pair, it was less than 2°C.

No clinical or non-clinical trials have been conducted to rule out the possibility of stent migration at an over 3.0 Tesla strength field. It is not recommended to perform MR tests with over 3.0 Tesla strength.

MR imaging quality may be affected if the area of interest is adjacent or very close to the stent place.

### 7. Drug interaction

Although specific clinical information is not available, drugs that act through the same binding protein (FKBP) as sirolimus may interfere with its efficacy. In turn, potent CYP3A4 inhibitors (such as ketoconazole) may cause an increase in exposure of sirolimus up to levels associated with systemic effects, especially when multiple stents are implanted. One should also keep in mind the systemic presence of sirolimus when treating a patient concomitantly with immunosuppression therapy.

### 8. Possible Adverse Effects / Complications

Possible associated adverse effects and/or complications that may arise before, during or after the procedure include the following:

[1] These data were obtained bibliographically, by searching for marketed stents with the same composition. Drug interaction

- Death
- Acute myocardial infarction
- Unstable angina
- Ictus / embolism / thrombosis
- Emergency bypass surgery
- Endocarditis
- Allergic reaction to the metals, to drugs, antiplatelet agents/anticoagulant agents or contrast media
- Artery spasm
- Hemodynamic deterioration
- Arrhythmia
- Aneurysm or pseudo-aneurysm
- Infections
- Total occlusion of the artery / thrombosis / stent embolization
- Reocclusion of the treated area: re-stenosis
- Perforation or dissection of the treated area
- Arteriovenous fistula
- Hypotension/hypertension
- Local haemorrhaging with haematoma at the insertion site

Possible adverse events associated with daily administration of sirolimus not previously indicated are:

- Anaemia
- Arthralgia
- Diarrhoea
- Interstitial lung disease
- Abnormal liver function
- Hypercholesterolaemia
- Hypertriglyceridaemia
- Hypersensitivity, including anaphylactic reactions
- Hypopotaemia
- Leucopenia
- Lymphoma and other malignant neoplasms
- Thrombocytopaenia

## 9. Individualization of treatment

Before using the product, the risks and benefits should be considered for each patient. Factors that should be kept in mind include an evaluation of the risk of using antiplatelet therapy. Special attention should be paid to patients with gastritis or recently active peptic ulcers.

## 10. How to Use

### 10.1 Required equipment

- Heparinised normal saline solution.
- Contrast medium (use contrast media indicated for intravascular use).
- Proper medical prescription to the procedure.
- 0.014" guidewire; do not use any other size of guidewire.
- Introducer with haemostatic valve whose size is equal or larger than 5F (except for 4 and 4.5 mm stent diameters, where 6F is recommended). Do not use smaller sizes as these might damage or leave useless the catheter.
- Guiding catheter or long introducer.
- Haemostatic valve.
- Three-way stopcock.
- Extension lines
- Inflation device (manual pump with built-in manometer).
- Several standard 10-20cc syringes with saline solution for washing the system.

### 10.2 System Preparation

- The stent diameter must coincide with the reference vessel diameter, as its length with lesion length. It is preferably to be a slightly longer to ensure full lesion coating.
  - Remove the catheter from the protective dispenser. Check that it is the right size. Remove the protective pod and stylet from the balloon. Wash the guide lumen with sterile saline solution by means of injection with a syringe from the tip until it outflows from the guidewire insertion port. Do not dry with gauze.
  - Inspect the stent to ensure it is located between the radiopaque markers; if it is noted that it has moved, do not use it.
  - Attach a three-way stopcock to an extension line and the extension line to the luer hub of the delivery system. Optionally, the inflation pump can be connected to the three-way stopcock.
  - Close the air passage through the balloon.
  - Attach a 10-20ml syringe with a third part of saline solution to the three-way stopcock.
  - Open the three-way stopcock connection between the syringe and the catheter.
  - Withdraw the syringe plunger upwards allowing the air bubbles exiting to the liquid (it is recommended that the syringe be in a vertical position).
  - When bubbles stop coming in, leave a meniscus of solution in the key and close the three-way stopcock at the catheter part. Remove the syringe.
  - Connect the inflation device with 1/3 of 1:1 contrast medium and saline solution mixture (duly purged), to the three-way stopcock and the extension line. Do not apply pressure as this might dislodge the stent.
- DANGER: If bubbles do not stop coming from the catheter into the syringe after 1 minute of negative pressure, this may be a clear indication that the balloon catheter has leaks, is broken or the connections between the syringe and the three-way stopcock are not properly sealed. If bubbles continue to be seen after checking the connections, do not use the device. Return the device to the manufacturer or distributor for its inspection.

### 10.3 Implantation Technique

#### Preliminary steps

- Insert the introducer according to the manufacturer's instructions.
- Connect the haemostatic valve to the guiding catheter hub or long introducer, maintain the valve closed. Insert the guiding catheter through the introducer or directly through the long introducer according to the manufacturer's instructions and make sure it is located at the peripheral artery entrance. Inject contrast liquid through the guiding catheter or the long introducer to ensure it is correctly positioned giving access to the artery.
- Open the haemostatic valve for 0.014" guidewire passage through. Once the guide crosses the key valve, seal it to avoid blood loss. Place the 0.014" guidewire through the lesion according to percutaneous intervention techniques, by means of fluoroscopy techniques to determine its position all times. Use a contrast liquid infusion to determine that the guidewire is properly positioned.

#### Implantation technique

- Pre-dilatation with a peripheral dilatation balloon catheter is recommended. Follow the instructions of the balloon catheter manufacturer for lesion pre-dilatation. Then, proceed to stent implantation.
- Insert the proximal end of the guidewire into the distal tip of the delivery system.
- Open the haemostatic valve again to allow stent passage through the valve. Once the stent has been run through it, close the valve. If any resistance is noted, do not advance the delivery system through the haemostatic key. Avoid the haemostatic key valve damaging the catheter shaft; which may subsequently affect the balloon in-deflation.
- Carefully advance the stent system on the guidewire through the guiding catheter until it reaches the lesion to treat, by means of fluoroscopic imaging. Ensure that the guidewire comes out of the balloon catheter guide insertion port (approximately 25cm from the tip of the balloon catheter).
- Place the stent across the lesion situating the radiopaque markers on both sides of the lesion. Use high-resolution fluoroscopy to verify that the stent has not been damaged or shifted during its advance before its expansion.

- Apply pressure for 15-30 seconds depending on the lesion, to completely expand the stent. Do not exceed the rated burst pressure (RBP) indicated at the label and the compliance curve.
- For an optimal expansion, the stent must be in contact with the artery, otherwise the stent could move when the balloon is removed. The balloon diameter must be slightly higher than the vessel diameter to allow stent elastic recoil once the balloon is deflated. Consult the compliance table for proper expansion, due to too high size may also induce artery dissection.
- Withdraw the plunger of the inflation device to deflate the balloon. Keep up the negative pressure for between 15 and 30 seconds depending on the balloon size. Ensure that the balloon is fully deflated (by means of fluoroscopy) before moving the catheter.
- With negative pressure in the inflation device and with the guidewire still placed, withdraw the catheter from the lesion. Keep the guide across the dilated stenosis.
- Carry out an angiography through the guiding catheter to confirm suitable stenting.
- Keep the guidewire through the dilated stenosis for 30 minutes after angioplasty. In case of any thrombus observed around the stent, administer thrombolytic agents. When the angiography confirms dilatation, carefully withdraw the guidewire and the balloon catheter through the connector.
- Withdraw the guiding catheter through the introducer.
- Leave the introducer placed until the haemodynamic profile returns to normality. Perform suture as usually.

#### Further dilatation of stented segments

- If the stent size results inadequate with respect to vessel diameter, stent post-dilatation technique may be applied with a larger balloon. To do this, a balloon catheter is inserted into the stented area (it is not needed to be compliant) and dilatation is performed according to the manufacturer's recommendations.
- As a further recommendation, a "small stent" must not be post-dilated over 3 mm, a "medium stent" must not be post-dilated over 4.25mm and a "large stent" must not be post-dilated over 5.25mm.

#### Unexpanded stent withdrawal

- If the stent does not cross the lesion, it may need to be removed without expanding. If any resistance is not detected, check that the guiding catheter is coaxially situated with respect to the stent and carefully withdraw the stent towards the guiding catheter. If any unusual resistance is noticed during its withdrawal, extract the whole system as a single unit by means of fluoroscopic imaging according to the indicated instructions in part 5.3 of this leaflet.

### 10.4 Antiplatelet therapy

- The optimal duration for double antiplatelet treatment is unknown, but clinical practice guidelines recommend administering aspirin indefinitely with clopidogrel for at least 6 months, stretching the treatment with clopidogrel up to 12 months in patients with low risk of bleeding.
- It is very important that the patient follows recommendations related to antiplatelet treatment after the intervention. Early interruption of the antiplatelet medication described could increase the risk of thrombosis, myocardial infarction or death. Before carrying out a percutaneous peripheral intervention (PPI), if a surgical or dental operation is planned for which it is necessary to interrupt antiplatelet treatment, the doctor and the patient should consider whether a drug-eluting stent and the associated antiplatelet treatment is the appropriate choice for the PPI. After the PPI, if carrying out a surgical or dental operation were to be prescribed, the risks and benefits of said operation should be weighed against the possible risk of an early suspension of antiplatelet treatment.

### 11. Warranty

The device and all its components have been designed, manufactured, tested and packed with a maximum level of precaution. LVD Biotech warrants the device until its expiration date, provided the packaging is not broken, tampered with or damaged.

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© iVascular - cod. MP45020 Rev. 0 may 2016