# essential

EN	INSTRUCTIONS FOR USE. PACLITAXEL ELUTING CORONARY BALLOON DILATATION CATHETER .....	3
ES	INSTRUCCIONES DE USO. CATÉTER BALÓN DE DILATACIÓN CORONARIO CON LIBERACIÓN DE PACLITAXEL .....	9
PT	INSTRUÇÕES DE USO. CATETER DE DILATAÇÃO CORONÁRIO COM LIBERAÇÃO DE PACLITAXEL .....	15
IT	ISTRUZIONI PER L' USO. CATETERE A PALLOCINO DI DILATAZIONE CORONARIA CON RILASCIO DI PACLITAXEL ....	21
FR	INSTRUCTIONS D'UTILISATION. CATHETER BALLONNET DE DILATATION CORONAIRE AVEC LIBERATION DE PACLIITAXEL	27
DE	GEBRAUCHSHINWEISE. PACLITAXEL-FREISETZENDER KORONARER DILATATIONSBALLONKATETHER.....	33
RU	ИНСТРУКЦИЯ ПО ПРИМЕНЕНИЮ. БАЛЛОННЫЙ КАТЕТЕР ДЛЯ КОРОНАРНОЙ ДИЛАТАЦИИ С ВЫСВОБОЖДЕНИЕМ ПАКЛИТАКСЕЛА .....	39
CS	NÁVOD K POUŽITÍ. BALÓNKOVÝ DILATAČNÍ KORONÁRNÍ KATÉTR S UVOLŇOVÁNÍM PACLITAXELU .....	45
EL	ΟΔΗΓΙΕΣ ΧΡΗΣΗΣ. ΚΑΘΗΤΗΡΑΣ ΜΕ ΜΠΑΛΟΝΙ ΕΚΛΟΥΣΗΣ ΠΑΚΛΙΤΑΞΕΛΗΣ ΣΤΕΦΑΝΙΑΙΩΝ ΑΡΤΗΡΙΩΝ .....	51
PL	INSTRUKCJA OBSŁUGI. CEWNIK BALONOWY UWALNIAJĄCY PAKLITAKSEL .....	57

## INSTRUCTIONS FOR USE

### PACLITAXEL ELUTING CORONARY BALLOON DILATATION CATHETER

#### 1. Description of the product

The paclitaxel eluting coronary balloon dilatation catheter is a rapid exchange catheter, also known as RX because it has a proximal simple lumen configuration and a distal coaxial dual lumen, intended for coronary arteries percutaneous transluminal angioplasties.

At the distal end of the catheter, before the tip, is the balloon or segment inflatable to different pressures which dilates the artery when inflated by means of an infusion of contrast liquid inside it. The balloon is designed to provide an inflatable segment of known diameter and length at the pressures recommended. The catheter ends in a cone-shaped tip.

The body of the catheter features a combination of a single lumen at the proximal end and dual lumen at the distal end.

- The single lumen at the proximal end allows contrast liquid to pass through between the proximal luer-lock connector and the balloon in order to inflate the balloon. This lumen is connected to the distal contrast liquid lumen at the end of the balloon.
- The coaxial dual lumen at the distal end has two functions:
  - a lumen to carry the contrast liquid used to inflate the balloon, connected to the proximal lumen
  - the guide wire passes through the other lumen to facilitate and allow the passage of the catheter towards and through the stenosis to be dilated.

There are marks in the catheter's shaft that help the operator to estimate the catheter's position in relation to the tip of the guide catheter (the marker located closer to the catheter connector is for femoral guide catheters and the farthest one is for brachial guide catheters)

The maximum diameter of the guide wire must be no more than 0.36 mm = 0.014".

The balloon is coated with the drug Paclitaxel, included in a physiologically innocuous matrix. The drug's dose is 3 micrograms / mm<sup>2</sup> of balloon surface and it is intended to avoid cellular proliferation, consequently decreasing the re-intervention rate.

The drug is released from the balloon by means of a rapid inflation in a way so that a high dose of the drug is released in a very short period of time. For releasing the drug onto the artery wall, the inflation process must last from 30 seconds to 1 minute, but dilation of the lesion must be optimized by using longer inflation times at the operator's discretion.

The catheter's useful length is 142 cm.

The product has the following balloon diameters and lengths.

		BALLOON LENGTH (mm)					
		10	15	20	25	30	40
BALLOON DIAMETERS (mm)	1,5	X	X	X		X	
	2,0	X	X	X	X	X	X
	2,25	X	X	X	X	X	X
	2,5	X	X	X	X	X	X
	2,75	X	X	X	X	X	X
	3,0	X	X	X	X	X	X
	3,25	X	X	X	X	X	X
	3,5	X	X	X	X	X	X
	3,75	X	X	X	X	X	X
	4,0	X	X	X	X	X	X
	4,25	X	X	X	X	X	X
4,5	X	X	X	X	X	X	

The balloon is preset to reach different diameters at different pressures (according to the compliance curve shown on the primary packaging).

To inflate the balloon, the proximal luer-lock connector must be connected to a manual pump fitted with a

**INSTRUCTIONS FOR USE****PACLITAXEL ELUTING CORONARY BALLOON DILATATION CATHETER**

manometer. The balloon is inflated to predictable diameters with the pressure shown by the manometer. Two radiopaque markers are fitted at each end of the balloon in order to mark its length and help the user to ascertain its whereabouts while inside the patient.

At the distal end of the balloon is the tip, which is rounded and atraumatic in shape in order to avoid damaging the arteries while it is being advanced.

**Contents**

- One rapid-exchange, Paclitaxel-releasing balloon catheter, covered by a protector over the balloon and with a guide lumen protection stylet. All the set is introduced in a dispenser to avoid damaging the catheter, and it is packed into a sterile bag.
- One card with the compliance curve showing the working range of pressures
- One leaflet with instructions for use

**2. Indications**

The Essential paclitaxel-eluting balloon catheter is indicated for the dilatation of stenosis or coronary artery or bypass grafts occlusions, including small vessels, as well as for residual stenosis after treatment with balloon or endoprosthesis and pre- and post-dilation of coronary endovascular prosthesis, in order to improve myocardial perfusion.

**3. Contraindications**

- Treatment of the left stem (first section of the left coronary artery)
- Coronary artery spasm without significant stenosis
- Patients showing angiographic evidence of the existence of thrombus.
- It must not be used in pregnant women or during lactation or in patients with known hypersensitivity to paclitaxel.

**4. Warnings**

- The device must be used by physicians with experience and trained in the percutaneous transluminal coronary angioplasty technique (PTCA).
- The coronary percutaneous intervention must be performed in hospitals equipped with the emergency facilities to perform open surgery operations or in hospitals with easy access to other hospitals where this surgery may be performed.
- This is a single-use product. It must not be resterilized or reused once the procedure is finished. Reusing this product with another patient may cause cross contamination, infection or transmission of infectious diseases from one patient to the next. Reusing the product may cause alterations thereof and limit its effectiveness.
- The product is supplied sterile. Check the expiry date and do not use products which are past this date.
- Inspect the pack before opening it. If there are any defects or the pack is damaged, throw the product away.
- Use aseptic techniques when the product is taken out of the primary packaging.
- Select the balloon diameter depending on the diameter of the artery to be treated. Do not insert a balloon with a diameter larger than that of the artery being treated.
- Do not touch the balloon directly or dry it with gauze, since doing so may cause the drug coating to come off.
- Do not expose the product to organic solvents such as alcohol, or to contrast media not indicated for intravascular use.
- Do not use air or any other gaseous media to inflate the balloon.
- Advance the product on the guide wire using fluoroscopy. Do not allow the product to advance without the guide wire inside it.
- Do not manipulate, advance or withdraw the catheter or the guidewire, when the balloon is inflated.
- In case a post-dilatation is necessary after use of product, it must be made with a conventional dila-

**INSTRUCTIONS FOR USE****PACLITAXEL ELUTING CORONARY BALLOON DILATATION CATHETER**

tation balloon, one without the drug, since the same vessel segment must not be dilated with more than one balloon with drug.

- Drug eluting stents must not be implanted in the same segment that has been previously treated with the balloon with drug, since it is possible for overdose or interaction among the active principles to occur.
- Published in vitro and in vivo studies prove the mutagenic effect of Paclitaxel. The carcinogenesis risk potential has not been studied. The published in vivo studies prove the teratogenic effect, therefore it is not recommended for pregnant women to use this product.

**5. Precautions**

- Administer suitable medical therapy to the patient: anticoagulants, vasodilators, etc., in accordance with the procedure for inserting coronary intravascular catheters.
- Before use, check the size and suitability of the catheter for its intended use.
- Use the product with precautions in complex lesions, since the abrasive nature of said lesions may alter the drug coating. In these cases it is recommended to predilate the lesion previously with a conventional, drugless balloon.
- Extreme care must be exercised when the device is taken off from the pack and passed through the hemostatic valve to make sure that the drug coating is not damaged.
- Proceed with extreme care so as not to damage the catheter while advancing it.
- The balloon must be inflated with a mixture of saline solution and contrast liquid (preferably 50/50, though the proportion of contrast liquid can be reduced in large balloons).
- If any resistance is noted during advance, stop and determine its cause before continuing.
- Do not exceed the rated burst pressure (RBP) as the balloon could burst.
- If any resistance to withdrawing the catheter is noted, it is recommended that the balloon catheter, the guide wire and inserter be removed as a single unit.
- Store in a cool, dry place away from direct sunlight.
- After use, this product may represent a biohazard. Handle and dispose of it in accordance with accepted medical practices and pertinent local, state or federal laws and regulations.

**6. Possible Adverse Effects / Complications**

The possible adverse effects and/or complications which might arise before, during or after the procedure include the following:

- Death
- Acute myocardial infarction
- Unstable angina
- Ictus / embolism / thrombosis
- Haemodynamic deterioration
- Acute elastic recoil
- Arrhythmia
- Aneurysm or pseudoaneurysm
- Infections
- Total occlusion of the artery
- Reocclusion of the area treated: restenosis
- Spasm
- Perforation or dissection of the area treated
- Arteriovenous fistula
- Local haemorrhaging with haematoma at the access site
- Hypo/hypertension
- Allergic or immune reaction to the drug.
- Histological changes in the arterial wall.

**INSTRUCTIONS FOR USE****PACLITAXEL ELUTING CORONARY BALLOON DILATATION CATHETER**

---

**7. Pharmacological interaction**

The potential interactions of paclitaxel with the drugs administered concomitantly have not been investigated. Drug interactions at the systemic levels of paclitaxel are detailed in the product data sheet but, considering that the amount of paclitaxel on the balloon surface corresponds approximately to a few hundredths of the amount normally used in the antineoplastic treatment, it is considered very unlikely for interactions with other drugs to occur.

Nevertheless, caution must be exercised when CYP3A4 or CYP2C8 substrates (such as terfenadine, cyclosporin, lovastatin, midazolam y ondansetron) or drugs with high PPB (Plasma Protein Binding) are administered concomitantly.

**8. How to Use****8.1 Equipment Required**

- Heparinised normal saline solution
- Contrast medium (use contrast media suitable for intravascular use)
- Medication suitable for the procedure
- 0.014" guide wire; do not use any other size of guide wire
- Inserter with a haemostatic valve
- Guide catheter
- Haemostasis valve
- Three-way stopcock
- Line extensions
- Inflating device (manual pump with built-in manometer)
- Several standard 10-20cc syringes with saline solution for washing the system

**8.2 Preparing the Catheter**

- In case of occlusion or very calcified lesions, it is recommended to predilate the lesion with a conventional balloon, without drug.
- Remove the catheter from the protective dispenser. Check that it is the right size. Gently remove the protective pod and stylet from the balloon.
- Purge the system of air.
  - Couple a three-way stopcock to the luer connector on the balloon catheter (optionally, to improve the mobility of the system, an extension can be coupled to the balloon connector and then the three-way stopcock).
  - Seal the air passage through the balloon.
  - Couple a 10-20ml syringe with a third part of saline solution to the three-way stopcock.
  - Open the connection in the three-way stopcock between the syringe and the catheter
  - With the syringe in a vertical position, withdraw the plunger upwards allowing the air bubbles to exit to the liquid.
  - When bubbles stop coming in, leave a meniscus of solution in the key and close the three-way stopcock at the catheter end. Remove the syringe.
  - Connect the inflating device with a 1/3 mixture of contrast liquid and saline solution (duly purged) to the three-way stopcock.

**DANGER:** If bubbles do not stop coming from the catheter into the syringe after 1 minute of negative pressure, this may be a clear indication that the balloon catheter has leaks, is broken or the connections between the syringe and the three-way stopcock are not properly sealed. If bubbles continue to be seen after checking the connections, do not use the device. Return it to the manufacturer or distributor for inspection.

**INSTRUCTIONS FOR USE****PACLITAXEL ELUTING CORONARY BALLOON DILATATION CATHETER****8.3 Technique for Insertion / Treatment**

- Insert the introducer according to the manufacturer's instructions.
- Connect the haemostasis valve to the guiding catheter connector, keeping the valve closed. Insert the guide catheter through the introducer, according to the manufacturer's instructions and making sure it is located at the entrance to the coronary artery. Inject contrast liquid through the catheter to make sure it is correctly positioned to give access to the coronary artery.
- Open the haemostasis valve to run the 0.014" guide wire through. Once the guide wire is through the valve, seal it to avoid loss of blood. Pass the 0.014" guide wire through the lesion according to percutaneous intervention techniques, with the aid of fluoroscopy to determine its position at all times. Use contrast liquid infusion to ensure that the guide wire is properly positioned.
- Insert the proximal end of the guide wire into the distal end of the balloon catheter.
- Open the hemostatic valve to allow passage of the balloon tip through the valve. Once the balloon is through, turn off the valve. If resistance is met, do not advance the dilation catheter through the hemostatic valve. Prevent the hemostatic valve from damaging the catheter body, which may later influence the inflation or deflation of the balloon. Take special care when passing the balloon catheter through the hemostatic valve, so as not to alter the drug coating.
- Advance carefully the balloon catheter along the guidewire, through the guide catheter and with the help of fluoroscopy, until the area to be treated is reached. Make sure that the guidewire goes out the insertion port of the balloon catheter guide (approximately 25 cm from the balloon catheter tip). In case the balloon catheter chosen cannot be placed inside the lesion, choose a smaller size of a drugless conventional balloon in order to predilate the lesion.
- Locate the balloon in the desired position for dilation.
- Apply pressure to the inflation device so that the balloon inflates. Do not exceed the rated burst pressure (RBP) on the label and the compliance curve.
- Keep the pressure for 30-60 seconds for an optimal release of the drug. Most of the drug is released in the first 30 seconds of inflation but, in order to optimize the dilation of the lesion, it is possible to use longer inflation times greater than 1 minute at the operator's discretion.
- Draw the plunger on the inflation device back to deflate the balloon. Keep up the negative pressure for between 15 and 30 seconds depending on the size of the balloon. Make sure the balloon is fully deflated (by means of fluoroscopy) before moving the catheter.
- With negative pressure in the inflation device and with the guide wire in position, remove the catheter carefully until it is pulled out through the introducer. Keep the guide wire in place across the dilated stenosis.
- Carry out an angiography through the guiding catheter to confirm dilation.
- Keep the guide wire through the dilated stenosis for 10 minutes after angioplasty. When the angiography confirms dilation, carefully withdraw the guide wire.
- Withdraw the guiding catheter through the introducer
- Leave the introducer in place until the haemodynamic profile returns to normal. Stitch in the usual way.
- Double antiplatelet therapy must be administered for a minimum of four weeks after the intervention, in the case of lesions without stent. In lesions with stent, a double antiplatelet treatment for 12 weeks is recommended. However, the treatment may be prolonged at the operator's discretion.

**Procedure for exchanging the dilation catheter**

In case it is necessary to change the balloon, the rapid exchange dilation catheter Essential is designed to allow one rapid exchange by only one operator. To make the exchange:

- Loosen the haemostatic valve.
- Hold the guide wire and the haemostatic valve in one hand and the balloon in the other.
- Keep the guide wire in position in the coronary artery; with the guide wire in a fixed position, withdraw the balloon catheter from the guide catheter.

**INSTRUCTIONS FOR USE**  
**PACLITAXEL ELUTING CORONARY BALLOON DILATATION CATHETER**

- Withdraw the catheter until it reaches the opening in the guiding catheter lumen (about 25cm from the end of the balloon). Carefully withdraw the distal portion of the balloon catheter, keeping the guide wire in place in the lesion. Close the haemostatic valve.
- Prepare the new balloon catheter to be used as described above.
- Fit the new balloon on the guide wire. Make sure the proximal portion of the guide wire exits through the balloon catheter guide insertion port to about 25cm from the distal end of the balloon.
- Open the haemostatic valve and advance the balloon while keeping the guide wire in the same position in the coronary artery. Care must be taken not to bend or twist the balloon catheter around the guide wire.
- Advance the balloon catheter to the end of the guide catheter. Continue with PCTA using the method described above.

**9. Warranty**

The device and all its components have been designed, manufactured, tested and packed with a maximum level of precaution. LVD Biotech warrants the device until its expiration date, provided the packaging is not broken, tampered with or damaged.

**STANDARDIZED SYMBOLS FOR USE IN LABELLING MEDICAL DEVICES**  
**(EN 980:2008 STANDARD)**

Catalogue reference

Batch no.

Expiry date

Keep away from sunlight

Manufacturer

Sterilized with ethylene oxide

Limitation in storage temperature

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