# luminor 14 m

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**INSTRUCTIONS FOR USE. PACLITAXEL ELUTING  
PERIPHERAL BALLOON DILATATION CATHETER, FOR 0.014" GUIDEWIRE**

**1. Description of the Product**

The Paclitaxel eluting dilatation balloon is a coaxial catheter (double lumen from the connector to the tip, also called OTW), designed for percutaneous transluminal angioplasties of small peripheral arteries.

The catheter has a coaxial dual-lumen body from the connector to the tip.

On the catheter's distal part, before the tip, there is the balloon with the drug or segment inflatable at different pressures that dilates the artery upon inflating by means of infusion of contrast fluid inside it, the connector is Y-shaped and it has two entry ports.

- The side port allows contrast medium passage to dilate the balloon
- The straight port is the port for guidewire passage. The guidewire is introduced in the connector's entry port and it goes out the catheter's tip to take it to the lesion.

Moreover, the lumen inner of the catheter is the one through which the guide wire runs, while the outer lumen is intended for contrast medium passage.

The maximum diameter of the guide wire must not exceed 0.36mm = 0.014 inches.

The balloon is coated with the drug Paclitaxel, included in a physiologically innocuous matrix. The drug's dose is 3 micrograms / mm<sup>2</sup> of balloon surface and it is intended to avoid cellular proliferation, consequently decreasing the re-intervention rate.

The drug is released from the balloon by means of a rapid inflation in a way so that a high dose of the drug is released in a very short period of time. For releasing the drug onto the artery wall, the inflation process must last from 30 seconds to 1 minute, but dilation of the lesion must be optimized by using longer inflation times at the operator's discretion.

The useful length of the catheter is 100cm or 150cm depending on the item reference.

The product has the following balloon diameters and lengths:

		BALLOON DIAMETERS (mm)					
		1,5	2,0	2,5	3,0	3,5	4,0
BALLOON LENGTH (mm)	40	X	X	X	X	X	X
	60	X	X	X	X	X	X
	80	X	X	X	X	X	X
	120		X	X	X	X	X
	150		X	X	X	X	
	200		X	X	X	X	

The balloon is intended to reach different diameters at different pressures (according to the compliance curve included on the primary packaging).

To inflate the balloon, the connector inflation port must be connected to an inflation device. The balloon is inflated to predictable diameters with the pressure shown by the manometer. At each end of the balloon there is a Pt/Ir radiopaque marker in order to mark its length and help the user to see the catheter while inside the patient.

On the distal end, there is the tip, it is rounded and atraumatic to avoid damaging the arteries during advancement.

**Contents**

- A paclitaxel-releasing dilatation balloon catheter OTW, covered by a protector over the balloon and a guide lumen protection stylet. All the set is introduced in a dispenser to avoid damaging the catheter, and it is packed into a sterile bag.
- One card with the compliance curve showing the nominal inflation pressure and the recommended maximum pressure
- One leaflet with instructions for use

**INSTRUCTIONS FOR USE. PACLITAXEL ELUTING  
PERIPHERAL BALLOON DILATATION CATHETER, FOR 0.014" GUIDEWIRE****2. Indications**

The device is indicated to dilate stenosis in femoral, popliteal and infra-popliteal arteries, with a reference diameter balloon from 1.5 and 4.0 mm and lengths from 40 up to 200 mm.

**3. Contraindications**

Inability to cross lesion with guide-wire

It must not be used in pregnant women or during lactation or in patients with known hypersensitivity to paclitaxel.

**4. Warnings**

- The device should only be used by experienced doctors familiar in the percutaneous transluminal angioplasty (PTA) technique.
- The percutaneous procedure must be carried out in hospitals equipped with emergency facilities to conduct open heart surgery or hospitals with easy access to other hospitals where this type of operation can be performed.
- The product is for single use. It must not be re-sterilised or re-used once the procedure has been completed. Re-using the product in another patient may lead to cross contamination, infection or transmission of infectious diseases from one patient to another. Re-use of the product may cause alterations to it and affect its effectiveness.
- The product is supplied sterile. Check the expiry date and do not use products which are past this date.
- Inspect the packaging before opening it. If there are any defects or the packaging is damaged, do not use the product.
- Use aseptic techniques when the product is taken out of the primary packaging.
- Select the balloon diameter according to the diameter of the artery to be treated. Do not insert a balloon with a diameter larger than the artery to be treated.
- Do not touch the balloon directly or dry it with gauze, since doing so may cause the drug coating to come off.
- Do not expose the product to organic solvents such as alcohol, or to contrast media not indicated for intravascular use.
- Do not use air or gaseous media to inflate the balloon.
- Advance the product on the guide wire using fluoroscopy. Do not allow the product to advance without the guide wire inside it.
- Do not manipulate, advance or retract either the catheter or the wire guide while the balloon is inflated.
- In case a post-dilatation is necessary after use of product, it must be made with a conventional dilatation balloon, one without the drug, since the same vessel segment must not be dilated with more than one balloon with drug.
- Drug eluting stents must not be implanted in the same segment that has been previously treated with the balloon with drug, since it is possible for overdose or interaction among the active principles to occur.
- Published in vitro and in vivo studies prove the mutagenic effect of Paclitaxel. The carcinogenesis risk potential has not been studied. The published in vivo studies prove the teratogenic effect, therefore it is not recommended for pregnant women to use this product.

**5. Precautions**

- Administer suitable medical therapy to the patient: anticoagulants, vasodilators, etc., in accordance with the procedure for inserting intravascular catheters.
- Use with standard 0.014-inch guide wires and standard introducers. If the introducer length goes above 25cm or the braided type is used, it may be necessary to increase its size to avoid catheter friction.
- Before use, check the size and suitability of the catheter for its intended use.
- Use the product with precautions in complex lesions, since the abrasive nature of said lesions may alter the drug coating. In these cases it is recommended to predilate the lesion previously with a con-

**INSTRUCTIONS FOR USE. PACLITAXEL ELUTING  
PERIPHERAL BALLOON DILATATION CATHETER, FOR 0.014" GUIDEWIRE**

ventional, drugless balloon.

- Extreme care must be exercised when the device is taken off from the pack and passed through the hemostatic valve to make sure that the drug coating is not damaged.
- Proceed with extreme care so as not to damage the catheter while advancing it.
- The balloon must be inflated with a mixture of saline solution and contrast medium (preferably 50/50, though the proportion of contrast medium can be reduced in large balloons, like 4mm diameter or longer ones).
- If any resistance is noticed during advance, stop and determine its cause before continuing.
- Do not exceed the rated burst pressure (RBP) as the balloon could burst.
- If any resistance to withdrawing the catheter is noted, it is recommended that the balloon catheter, the guide wire and introducer be removed as a single unit.
- Store in a cool, dry place away from direct sunlight.
- After use, this product may represent a biohazard. Handle and dispose it in accordance with accepted medical practices and pertinent local, state or federal laws and regulations.

**6. Possible Adverse Effects / Complications**

The possible adverse effects and/or complications which might arise before, during or after the procedure include the following:

- Death
- Ictus / embolism / thrombosis
- Haemodynamic deterioration
- Acute elastic recoil
- Arrhythmia
- Aneurysm or pseudo-aneurysm
- Infections
- Total occlusion of the artery
- Reocclusion of the treated area: restenosis
- Spasm
- Perforation or dissection of the treated area
- Arteriovenous fistula
- Local bleeding with haematoma at the access site
- Hypo/hypertension
- Allergic or immune reaction to the drug.
- Histological changes in the arterial wall.

**7. Pharmacological interaction**

The potential interactions of paclitaxel with the drugs administered concomitantly have not been investigated. Drug interactions at the systemic levels of paclitaxel are detailed in the product data sheet but, considering that the amount of paclitaxel on the balloon surface corresponds approximately to a few hundredths of the amount normally used in the antineoplastic treatment, it is considered very unlikely for interactions with other drugs to occur.

Nevertheless, caution must be exercised when CYP3A4 or CYP2C8 substrates (such as terfenadine, cyclosporin, lovastatin, midazolam y ondansetron) or drugs with high PPB (Plasma Protein Binding) are administered concomitantly.

**8. How to Use****8.1 Equipment Required**

- Heparinized normal saline solution
- Contrast medium (use contrast media suitable for intravascular use)
- 0.014" guide wire; do not use any other size of guide wire
- Introducer with a hemostatic valve of the size specified on the label. Do not use smaller sizes as

## INSTRUCTIONS FOR USE. PACLITAXEL ELUTING PERIPHERAL BALLOON DILATATION CATHETER, FOR 0.014" GUIDEWIRE

these might damage the catheter.

- Three-way stopcock
- Line extensions
- Inflation device (manual pump with built-in manometer)
- Several standard 10-20cc syringes with saline solution for washing the system

### **8.2 Preparing the Catheter**

- In case of occlusion or very calcified lesions, it is recommended to predilate the lesion with a conventional balloon, without drug.
- Remove the catheter from the protective dispenser. Check that it is the right size. Gently remove the protective sheath and stylet from the balloon.
- Connect a 10ml syringe containing sterile saline solution to the port for the guide wire (straight part of the connector) and irrigate the lumen until the liquid comes out of the tip.
- Purge the whole system of air.
  - Attach a three-way stopcock to the contrast liquid port (the side part of the connector). An extension line could be optionally placed between the balloon connector and the three-way stopcock.
  - Close the air passage through the balloon
  - Attach a 10-20ml syringe with a third part of saline solution to the three-way stopcock.
  - Open the passage in the three-way stopcock between the syringe and the catheter
  - With the syringe in a vertical position, withdraw the plunger upwards allowing the air bubbles to exit to the liquid.
  - When bubbles stop coming in, close the three-way stopcock at the catheter end and remove the syringe

**DANGER:** If bubbles do not stop coming from the catheter into the syringe after 3 minutes of negative pressure, this may be a clear indication that the balloon catheter has leaks, is broken or the connections between the syringe and the three-way stopcock are not properly closed. If bubbles continue to be seen after checking the connections, do not use the device. Return it to the manufacturer or distributor for inspection.

### **8.3 Technique for Insertion / Treatment**

- Position the 0.014" guide wire across the lesion using PTA techniques, using fluoroscopy to determine its position at all times.
- Insert the proximal end of the guide wire into the distal tip of the balloon catheter.
- Advance carefully the balloon catheter along the guidewire, through the guide catheter and with the help of fluoroscopy, until the area to be treated is reached. Extreme care must be exercised when the balloon catheter is passed through the hemostatic valve to make sure that the drug coating is not damaged. In case the balloon catheter chosen cannot be placed inside the lesion, choose a smaller size of a drugless conventional balloon in order to predilate the lesion.
- Connect the inflation device (1/3 of its capacity of contrast liquid mixture: saline solution) to the three-way stopcock and open the passage between the balloon and the inflation device.
- Apply pressure to the inflation device so that the balloon inflates. Do not exceed the rated burst pressure (RBP) on the label and the compliance curve.
- Keep the pressure for 30-60 seconds for an optimal release of the drug. Most of the drug is released in the first 30 seconds of inflation but, in order to optimize the dilation of the lesion, it is possible to use longer inflation times greater than 1 minute at the operator's discretion.
- Draw the plunger on the inflation device back to deflate the balloon. Keep up the negative pressure for between 20 and 60 seconds depending on the size of the balloon. Make sure the balloon is fully deflated (by means of fluoroscopy) before moving the catheter.
- With negative pressure in the inflation device and with the guide wire in position, carefully retract the catheter until it can be removed through the introducer.
- Inject contrast liquid through the introducer to check the treated lesion before removing the guide wire.

**INSTRUCTIONS FOR USE. PACLITAXEL ELUTING  
PERIPHERAL BALLOON DILATATION CATHETER, FOR 0.014" GUIDEWIRE**

- Double antiplatelet therapy must be administered for a minimum of four weeks after the intervention, and the treatment may be prolonged at the operator's discretion.

**9. Warranty**

The device and all its components have been designed, manufactured, tested and packed with a maximum level of precaution. LVD Biotech guarantees the product up to its expiry date providing its packaging is not broken, tampered with or damaged.

**STANDARDIZED SYMBOLS FOR USE IN LABELLING MEDICAL DEVICES  
(EN 980:2008 STANDARD)**

Catalogue reference

Batch no.

Expiry date

Keep away from sunlight

Manufacturer

Sterilized with ethylene oxide

Limitation in storage temperature

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