
# Xperience

EN	INSTRUCTIONS FOR USE. BALLOON CATHETER FOR CORONARY PRE-DILATATION, WITH HYDROPHILIC COATING .....	2
ES	INSTRUCCIONES DE USO. CATÉTER BALÓN DE PREDILATACIÓN CORONARIO, CON RECUBRIMIENTO HIDROFÍLICO .....	7
PT	INSTRUÇÕES DE USO, CATETER BALÃO DE PRÉ-DILATAÇÃOCORONÁRIA, COM REVESTIMENTO HIDROFÍLICO .....	12
IT	ISTRUZIONI PER L'USO, CATETERE A PALLONCINO PER DILATAZIONE CORONARICA, CON RIVESTIMENTO IDROFILICO .....	17
FR	MODE D'EMPLOI. CATHÉTER À BALLONNET DE PRÉ-DILATATION CORONAIRE, AVEC REVÊTEMENT HYDROPHILE .....	22
DE	GEBRAUCHSHINWEISE. BALLONKATHETER FÜR KORONARE VORDILATATION, MIT HYDROPHILER BESCHICHTUNG .....	27
RU	ИНСТРУКЦИИ ПО ИСПОЛЬЗОВАНИЮ, БАЛЛОННЫЙ КАТЕТЕР С ГИДРОФИЛЬНЫМ ПОКРЫТИЕМ ДЛЯ ПРЕДИЛАТАЦИИ КОРОНАРНЫХ АРТЕРИЙ .....	32
CS	NÁVOD K POUŽITÍ. BALÓNKOVÝ KATÉTR S HYDROFILNÍM POVLAKEM PRO PREDILATACI KORONÁRNÍCH TEPEN .....	37
EL	ΟΔΗΓΙΕΣ ΧΡΗΣΗΣ. ΚΑΘΗΤΗΡΑΣ ΜΠΑΛΟΝΙ ΓΙΑ ΠΡΟΔΙΑΣΤΟΛΗ ΣΤΕΦΑΝΙΑΙΑΩΝ ΑΡΤΗΡΙΩΝ ΜΕ ΥΔΡΟΦΙΛΗ ΕΠΙΚΑΛΥΨΗ .....	42
PL	INSTRUKCJA UŻYTKOWANIA. CEWNIK BALONOWY Z POWŁOKĄ HYDROFILOWĄ DO ROZSZERZANIA NACZYŃ WIEŃCOWYCH .....	47

**INSTRUCTIONS FOR USE**

**BALLOON CATHETER FOR CORONARY PRE-DILATATION, WITH HYDROPHILIC COATING**

**1. Description of the Product**

The Xperience balloon catheter for coronary dilatation is a rapid exchange catheter, also referred to as RX, because it features a single proximal lumen and dual coaxial distal lumen, designed for percutaneous transluminal angioplasty of coronary arteries.

At the distal end of the catheter, before the tip, is the balloon or segment inflatable to different pressures which dilates the artery when inflated by means of an infusion of contrast liquid inside it. The balloon is designed to provide an inflatable segment of known diameter and length at the pressures recommended. The catheter ends in a cone-shaped tip.

The body of the catheter features a combination of a single lumen at the proximal end and dual lumen at the distal end.

- The single lumen at the proximal end allows contrast liquid to pass through between the proximal luer-lock connector and the balloon in order to inflate the balloon. This lumen is connected to the distal contrast liquid lumen at the end of the balloon.
- The coaxial dual lumen at the distal end has two functions:
  - a lumen to carry the contrast liquid used to inflate the balloon, connected to the proximal lumen
  - the guide wire passes through the other lumen to facilitate and allow the passage of the catheter towards and through the stenosis to be dilated.

There are marks on the body of the catheter to help calculate the position of the catheter in relation to the end of the guide catheter (the marker located nearer to the catheter adapter is for femoral guiding catheters and the one further away is for brachial guiding catheters).

The distal part of the catheter is covered in a durable hydrophilic coating to lubricate the catheter so that it can navigate through arteries.

The maximum diameter of the guide wire must be no more than 0.36 mm = 0.014".

The usable catheter is 142 cm long.

The balloon is preset to reach different diameters at different pressures (according to the compliance curve shown on the primary packaging).

To inflate the balloon, the proximal luer-lock connector must be connected to a manual pump fitted with a manometer. The balloon is inflated to predictable diameters with the pressure shown by the manometer. Two radiopaque markers are fitted at each end of the balloon in order to mark its length and help the user to ascertain its whereabouts while inside the patient.

At the distal end of the balloon is the tip, which is rounded and atraumatic in shape in order to avoid damaging the arteries while it is being advanced.

The XPERIENCE product is sold under different reference numbers / variants, depending on the dimensions of the balloon:

BALLOON DIMENSIONS		Lengths (mm)					
		10	15	20	25	30	40
Diameters (mm)	1.50	x	x	x		x	
	2.00	x	x	x	x	x	x
	2.25	x	x	x	x	x	x
	2.50	x	x	x	x	x	x
	2.75	x	x	x	x	x	x
	3.00	x	x	x	x	x	x
	3.25	x	x	x	x	x	x
	3.50	x	x	x	x	x	x
	3.75	x	x	x	x	x	x
	4.00	x	x	x	x	x	x
	4.25	x	x	x	x	x	x
	4.50	x	x	x	x	x	x

## INSTRUCTIONS FOR USE

### BALLOON CATHETER FOR CORONARY PRE-DILATATION, WITH HYDROPHILIC COATING

#### **Contents**

- One rapid-exchange balloon catheter covered by protection for the balloon and a stylet in the guide wire lumen, inserted into a dispenser and packed in a sterile bag
- One card with the compliance curve showing the working range of pressures
- One leaflet with instructions for use

#### **2. Indications**

- The Xperience balloon catheter is suitable for balloon dilatation of the stenotic portion of a coronary artery or of stenosis of bypass grafts replacing coronary arteries, in order to improve myocardial perfusion.

#### **3. Contraindications**

Unprotected left main coronary artery

Coronary artery spasm without significant stenosis

#### **4. Warnings**

- The product is single use. It must not be re-sterilised or re-used once the procedure has been completed. Re-using the product in another patient may lead to cross contamination, infection or transmission of infectious diseases from one patient to another. Re-use of the product may cause alterations to it and limit its effectiveness.
- The product is supplied sterile. Check the expiry date and do not use products which are past this date.
- Inspect the pack before opening it. If there are any defects or the pack is damaged, throw the product away.
- Use aseptic techniques when the product is taken out of the primary packaging.
- Select the balloon diameter depending on the diameter of the artery to be treated. Do not insert a balloon with a diameter larger than that of the artery being treated.
- Do not dry with gauze.
- Advance the product on the guide wire using fluoroscopy. Do not allow the product to advance without the guide wire inside it.
- Do not use air or gaseous media to inflate the balloon.
- Do not manipulate, advance or retract either the catheter or the wire guide while the balloon is inflated.
- Do not expose the product to organic solvents or to contrast media unsuitable for intravascular use.
- Percutaneous coronary intervention (PCI) must be carried out in hospitals equipped with emergency facilities to conduct open heart surgery or hospitals with easy access to other hospitals where this type of operation can be performed.

#### **5. Precautions**

- The device is for use by experienced doctors familiar with the percutaneous transluminal angioplasty (PTA) technique.
- Administer suitable medical therapy to the patient: anticoagulants, vasodilators, etc., in accordance with the procedure for inserting coronary intravascular catheters.
- Before use, check the size and suitability of the catheter for its intended use.
- Proceed with extreme care so as not to damage the catheter while advancing it.
- If any resistance is noted during advance, stop and determine its cause before continuing.
- If any resistance to withdrawing the catheter is noted, it is recommended that the balloon catheter, the guide wire and inserter be removed as a single unit.
- The balloon must be inflated with a mixture of saline solution and contrast liquid (preferably 50/50, though the proportion of contrast liquid can be reduced in large balloons).

**INSTRUCTIONS FOR USE****BALLOON CATHETER FOR CORONARY PRE-DILATATION, WITH HYDROPHILIC COATING**

- Do not exceed the rated burst pressure (RBP) as the balloon could burst.
- Store in a cool, dry place away from direct sunlight.
- After use, this product may represent a biohazard. Handle and dispose of it in accordance with accepted medical practices and pertinent local, state or federal laws and regulations.

**6. Possible Adverse Effects / Complications**

The possible adverse effects and/or complications which might arise before, during or after the procedure include the following:

- Death
- Acute myocardial infarction
- Unstable angina
- Ictus / embolism / thrombosis
- Coronary artery spasms
- Pseudoaneurysm
- Haemodynamic deterioration
- Arrhythmia
- Aneurysm or pseudoaneurysm
- Infections
- Total occlusion of the artery
- Reocclusion of the area treated: restenosis
- Perforation or dissection of the area treated
- Arteriovenous fistula
- Local haemorrhaging with haematoma at the access site
- Allergy to materials

**7. How to Use****7.1. Equipment Required**

- Heparinised normal saline solution
- Contrast medium (use contrast media suitable for intravascular use)
- Medication suitable for the procedure
- 0.014" guide wire; do not use any other size of guide wire
- Inserter with a haemostatic valve
- Guide catheter
- Haemostasis valve
- Three-way stopcock
- Line extensions
- Inflating device (manual pump with built-in manometer)
- Several standard 10-20cc syringes with saline solution for washing the system

**7.2. Preparing the Catheter**

- Remove the catheter from the protective dispenser. Check that it is the right size. Gently remove the protective pod and stylet from the balloon.
- Submerge the catheter in sterile saline solution to activate the hydrophilic coating.
- Purge the system of air.
- Couple a three-way stopcock to the luer connector on the balloon catheter (optionally, to improve the mobility of the system, an extension can be coupled to the balloon connector and then the three-way stopcock).
- Seal the air passage through the balloon.

**INSTRUCTIONS FOR USE****BALLOON CATHETER FOR CORONARY PRE-DILATATION, WITH HYDROPHILIC COATING**

- Couple a 10-20ml syringe with a third part of saline solution to the three-way stopcock.
- Open the connection in the three-way stopcock between the syringe and the catheter
- With the syringe in a vertical position, withdraw the plunger upwards allowing the air bubbles to exit to the liquid.
- When bubbles stop coming in, leave a meniscus of solution in the key and close the three-way stopcock at the catheter end. Remove the syringe.
- Connect the inflating device with a 1/3 mixture of contrast liquid and saline solution (duly purged) to the three-way stopcock.

**DANGER:** If bubbles do not stop coming from the catheter into the syringe after 1 minute of negative pressure, this may be a clear indication that the balloon catheter has leaks, is broken or the connections between the syringe and the three-way stopcock are not properly sealed. If bubbles continue to be seen after checking the connections, do not use the device. Return it to the manufacturer or distributor for inspection.

**7.3. Technique for Insertion / Treatment**

- Insert the introducer according to the manufacturer's instructions.
- Connect the haemostasis valve to the guiding catheter connector, keeping the valve closed. Insert the guide catheter through the introducer, according to the manufacturer's instructions and making sure it is located at the entrance to the coronary artery. Inject contrast liquid through the catheter to make sure it is correctly positioned to give access to the coronary artery.
- Open the haemostasis valve to run the 0.014" guide wire through. Once the guide wire is through the valve, seal it to avoid loss of blood. Pass the 0.014" guide wire through the lesion according to percutaneous intervention techniques, with the aid of fluoroscopy to determine its position at all times. Use contrast liquid infusion to ensure that the guide wire is properly positioned.
- Insert the proximal end of the guide wire into the distal end of the balloon catheter.
- Open the haemostasis valve again to let the end of the balloon pass through it. Once the balloon is through, close the valve. If resistance is met, do not advance the dilatation catheter through the haemostasis valve. Do not allow the haemostatic valve to damage the catheter body as this could then affect the inflation/deflation of the balloon.
- Carefully advance the balloon catheter on the guide wire, through the guide catheter and with the aid of fluoroscopy, until it reaches the area to be treated. Make sure the guide wire exits through the balloon catheter guide insertion port (about 25cm from the end of the balloon catheter). In the event that the balloon catheter selected cannot be fitted inside the lesion, choose a smaller size to predilate the lesion.
- Locate the balloon in the desired position for dilatation.
- Apply pressure to the inflation device so that the balloon inflates. Do not exceed the rated burst pressure (RBP) on the label and the compliance curve.
- Apply pressure for 15-30 seconds depending on the lesion.
- Draw the plunger on the inflation device back to deflate the balloon. Keep up the negative pressure for between 15 and 30 seconds depending on the size of the balloon. Make sure the balloon is fully deflated (by means of fluoroscopy) before moving the catheter.
- With negative pressure in the inflation device and with the guide wire in position, withdraw the catheter from the lesion. Keep the guide wire in place across the dilated stenosis.
- Carry out an angiography through the guiding catheter to confirm dilatation.
- Keep the guide wire through the dilated stenosis for 10 minutes after angioplasty. When the angiography confirms dilatation, carefully withdraw the guide wire and the balloon catheter through the connector.
- Withdraw the guiding catheter through the introducer

**INSTRUCTIONS FOR USE**

**BALLOON CATHETER FOR CORONARY PRE-DILATATION, WITH HYDROPHILIC COATING**

- Leave the introducer in place until the haemodynamic profile returns to normal. Stitch in the usual way.

**Procedure for exchanging the dilatation catheter**

- Should the balloon need changing, the Xperience rapid exchange dilatation catheter is designed to allow rapid exchange by a single operator. To exchange it proceed as follows:
  - Loosen the haemostatic valve.
  - Hold the guide wire and the haemostatic valve in one hand and the balloon in the other.
  - Keep the guide wire in position in the coronary artery; with the guide wire in a fixed position, withdraw the balloon catheter from the guide catheter.
  - Withdraw the catheter until it reaches the opening in the guiding catheter lumen (about 25cm from the end of the balloon). Carefully withdraw the distal portion of the balloon catheter, keeping the guide wire in place in the lesion. Close the haemostatic valve.
  - Prepare the new balloon catheter to be used as described above.
  - Fit the new balloon on the guide wire. Make sure the proximal portion of the guide wire exits through the balloon catheter guide insertion port to about 25cm from the distal end of the balloon.
  - Open the haemostatic valve and advance the balloon while keeping the guide wire in the same position in the coronary artery. Care must be taken not to bend or twist the balloon catheter around the guide wire.
  - Advance the balloon catheter to the end of the guide catheter. Continue with PCTA using the method described above.

**8. Warranty**

The product and all its components are designed, manufactured, tested and packed taking all reasonable precautions. LVD Biotech guarantees the product up to its expiry date providing its packaging is not broken, tampered or damaged.

STANDARDIZED SYMBOLS FOR USE IN LABELLING MEDICAL DEVICES  
(EN 980:2008 STANDARD)

			
Catalogue reference	Batch no.	Expiry date	Keep away from sunlight
			
Manufacturer	Sterilized with ethylene oxide	Limitation in storage temperature	Keep dry
			
Consult instructions for use	Do no re-use	Do not re-sterilize	Do not use if the pack is damaged

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