
# iVascular® iCover

EN	INSTRUCTIONS FOR USE. PERIPHERAL BALLOON EXPANDABLE PTFE COVERED STENT SYSTEM	03
ES	INSTRUCCIONES DE USO. SISTEMA DE STENT PERIFÉRICO CUBIERTO PTFE EXPANDIBLE POR BALÓN	05
PT	INSTRUÇÕES DE USO. SISTEMA DE STENT PERIFÉRICO REVESTIDO COM PTFE EXPANSÍVEL POR BALÃO	07
IT	ISTRUZIONI PER L'USO. SISTEMA DI STENT PERIFERICO RIVESTITO IN PTFE ESPANDIBILE CON PALLONCINO	09
FR	INSTRUCTIONS D'UTILISATION. SYSTÈME DE STENT PÉRIPHÉRIQUE RECOUVERT DE PTFE EXPANSIBLE PAR BALLONNE	11
DE	GEBRAUCHSHINWEISE. BALLONAUSSDEHNBARER PTFE-UMMANTELTER PERIPHERER STENTSYSTEM	13
RU	ИНСТРУКЦИИ ПО ИСПОЛЬЗОВАНИЮ. СИСТЕМА ПЕРИФЕРИЧЕСКОГО СТЕНТ С РАСШИРЯЮЩИМСЯ БАЛЛОННЫМ ПОКРЫТИЕМ ИЗ PTFE	16
CS	NÁVOD K POUŽITÍ. SYSTÉM BALÓNEK EXPANDOVATELNÝ PTFE POTAŽENÝ PERIFERNÍHO STENTEM	18
NL	GEBRUIKSAANWIJZING. PERIFEEER BALLON EXPANDEERBARE MET PTFE BEKLEDE STENTSYSTEEM	20

## ENGLISH

### 1. Product description

The covered peripheral stent system with balloon expansion iCover made by LVD Biotech SL is a stent manufactured using a Cobalt Chromium alloy called L605, which is subsequently covered internally and externally with ePTFE (expanded polytetrafluoroethylene) so that the stent is completely encapsulated. The covered stent system is presented preassembled in the positioning system that will expand it, by balloon, in the peripheral artery to be treated.

The stent is designed for different artery diameters by adapting its open cell design, alternating connection bridges. The stent is manufactured by laser cutting metal tubes. It is then subjected to various treatments that will give its surface a smooth and glossy finish. Once the stent has been manufactured, it is covered with ePTFE. The stent incorporates three radiopaque tantalum markers at each end, to indicate its position once implanted.

The stent release system consists of a double lumen coaxial balloon catheter from the proximal connector to the distal tip, also known as OTW, designed for percutaneous transluminal angioplasty of the peripheral arteries.

At its most distal end, the catheter has a balloon with a specific length and diameter, so it can be coupled to the different stent configurations. When the balloon is inflated by injecting contrast liquid into it, it dilates the stent and positions it against the artery. Then the balloon is deflated and withdrawn, and the stent remains permanently implanted.

The proximal connector is "Y-shaped" and has two entry ports:

- The lateral port for injecting contrast liquid to dilate the stent.
- The straight port is intended for guidewire insertion. The guidewire is inserted through the tip of the catheter and comes out by the connector port to lead the catheter to the lesion.

Similarly, there is a lumen for the guidewire and a lumen for the contrast liquid toward the balloon, in direct connection with each one of the ports of the luer or connector.

The diameter of the guidewire must not be greater than 0.89 mm = 0.353 inches.

The useful catheter length may be 80 cm or 140 cm, depending on the reference.

The stent can have the following diameters and lengths (pg. 23).

The balloon is predetermined to reach different diameters at different pressures (according to the compliance curve included in the primary packaging and pg 23).

To inflate the balloon, connect the inflation port of the connector to an inflation device, which at the same time is connected to a pressure gauge. With the pressure indicated by the pressure gauge, the balloon is inflated to predetermined diameters and dilates the covered stent. Each end of the balloon has a radiopaque marker to delimit the length and position of the stent, and to help the user know where it is as it advances inside the patient.

The tip is located at the distal end; it is rounded and atraumatic to avoid damaging the arteries during progress.

### 2. Contents

- A peripheral stent system composed of the covered stent preassembled on its release system (balloon catheter). The covered stent is protected by a protective sheath. The catheter is inside a dispenser positioned over a blister pack. The entire assembly is packaged in a sterile bag.
- A card with the compliance curve, indicating the range of working pressures.
- An implantation card.
- A leaflet with the instructions for use.

### 3. Intended use/Indications

Intended use:

- Restoring atherosclerotic lesions of the common and external iliac arteries and the renal arteries.
- Treatment of aneurysms and perforations, ruptures and acute fistulas.

Indications:

Patients presenting peripheral vascular system pathologies.

### 4. Contraindications

Patients with a clear diagnosis of allergy/hypersensitivity to heavy metals or ePTFE.

Patients with lesions that impede the complete inflation of the angioplasty balloon or the correct placement of the stent.

Presence of recent thrombus that has not dissolved, that is disorganised or of embolic material.

Unprotected vascular regions, exposed to the potential risk of external impact that could lead the stent to be crushed.

Aneurysms adjacent to the implantation site of the covered stent.

Patients that do not tolerate anticoagulant and/or antiplatelet treatment.

Threat of vital lateral branch occlusion.

### 5. Warnings

• Stents must only be implanted in hospitals with emergency facilities and an emergency team that can transfer the operation to open surgery, or, in their absence, hospitals with these services nearby.

• The product should only be used by physicians with experience in angioplasty techniques who are familiar with percutaneous intervention and stent implantation.

• This device is designed for single use in just one patient. DO NOT RESTERILISE OR REUSE. Reusing the product in another patient may cause cross-contamination, infection or transmission of infectious diseases. Reusing the product in another patient can alter it and compromise its efficacy.

• The product is supplied sterile. Check the expiry date and do not use products that have passed this date.

• Do not use the product if the packaging is open or damaged.

• Administer the appropriate medical treatment to the patient: anticoagulants, vasodilator, etc., according to the intravascular catheter insertion protocol.

• Use aseptic techniques when removing the product from its primary packaging.

• Do not dry with gauze.

• Do not expose the delivery device to organic solvents.

• **Before its implantation, moisten the whole catheter well, especially the stent, in saline solution for 20 seconds.**

• Do not use oily or high-viscosity contrast media that are not indicated for intravascular use.

• Do not use air or gaseous media to inflate the balloon; it must be inflated with a mixture of saline solution and contrast liquid (preferably 1:1 or 2:1 for large balloons).

• Select the appropriate product size in terms of diameter and length, according to the size of the lesion observed under fluoroscopy.

• Use materials that are compatible in accordance with the indication on the label, so as not to risk the integrity of the stent.

• Advance the product on the guidewire using fluoroscopic imaging. Do not allow the product to advance without the guidewire inside it.

• Insufficient expansion of the stent can cause it to migrate.

• For the maximum recommended postdilatation diameter, consult the table below. A larger postdilatation diameter than recommended could cause the fracture of the stent or the cover (pg. 23).

• Do not expose the covered stent to temperatures over 260°C, as the PTFE breaks down and can emit highly toxic decomposition products.

• Applying lasers to the surface of a stent covered with ePTFE, or around it, can damage the covered stent and generate toxic fumes.

• Do not superimpose stents of different compositions. Corrosion due to friction in superimposed prostheses has not been demonstrated. This device has not been tested in conditions of superimposition with stents or covered stents.

• Do not expose the product to hospital waste material.

### 6. Precautions

#### 6.1 Precautions during storage

- Store in a cool, dry place, away from direct sunlight.

#### 6.2 Precautions during preparation of the product

- Inspect the packaging. If it is damaged, do not use the product.
- The product label indicates the minimum size of the introducer with which the system must be used. Do not use introducers that are smaller than indicated. Moreover, if using a long or braided introducer, it may be necessary to increase one French size to ensure the success of the procedure.
- Remove the sheath and inspect the product before use to check it is in good condition, not curved or bent, and positioned between the radiopaque markers on the balloon. Dispose of any devices you are unsure of.
- Do not prepare or pre-inflate the release system before starting the procedure or extend the stent at the lesion site. Couple the inflation device with the diluted contrast medium in the inflation port.
- Purge the system and **moisten the stent in particular with saline solution for 20 seconds before implanting it.**
- Do not handle the stent with your fingers, as this can make the stent come out of the balloon.
- Do not try to separate the stent from its release system, if you notice that the stent moves or is not correctly positioned, do not use it.
- Before use, check the size and adapt the diameter of the stent to the diameter of the artery to dilate.

#### 6.3 Precautions during advancement - positioning - deployment - withdrawal of the product

- When treating multiple lesions, start with the most distal lesions, followed by the more proximal.
- Anatomical variations can complicate the procedure; precaution is recommended when advancing the endovascular system in tortuous and difficult configurations.
- Always advance the device under fluoroscopic visualisation. If you note resistance, do not continue until you have corrected the cause. If the stent cannot pass or reach the lesion, remove the whole system as an assembly, including the introducer.
- Before extending the stent, ensure it is in the correct location. Tense the release balloon catheter to prevent the stent from being positioned incorrectly. Do not inflate the balloon until you are sure it is in the correct position. Verify its position using fluoroscopy.
- Once you have started to expand the stent, do not attempt to withdraw it or change its position. Do not manipulate, advance or remove the catheter or the guidewire while the balloon is inflated. Only advance and withdraw the system if the balloon is fully deflated.
- Do not exceed the maximum recommended pressure indicated on the label and in the attached compliance curve.
- If the balloon bursts when you have started to expand the covered stent, take extreme care when removing it. Try to remove it together with the introducer as one single assembly. If it is not possible to remove the stent, leave the guidewire in position so that you can very carefully insert another semi-compliant balloon with the same diameter, to position the covered stent over the artery, while withdrawing the balloon and introducer as one single assembly. Once the stent is secure, deflate and remove the balloon and, if necessary, postdilate it again with another balloon. If it is possible to remove the stent, do so together with the balloon, guidewire and introducer as one single assembly.
- If you sense any kind of resistance during access to the lesion, when extracting the release system, or when withdrawing a stent that could not be implanted, the whole assembly must be withdrawn as a single unit:
  - Do not withdraw the delivery system inside the introducer since the stent may be released into the vasculature.
  - Position the proximal marker of the balloon distally to tip of the introducer.
  - Advance the guidewire forward in the anatomy, as far as possible.
  - Withdraw it until the introducer and the balloon catheter with the stent have been extracted together, as a single unit.
  - Lastly, withdraw the guidewire or, if you wish to start implantation again, leave the guidewire in place and position the other devices.
- Positioning the stent in a lateral branch or bifurcation can compromise future diagnostic or therapeutic procedures.
- Immediately after an implant, and if it is necessary to re-cross the stent, take great care to avoid pushing the stent with the guidewires or balloons to avoid dislodging it.

- After use, this product may pose a biological hazard. Handle and dispose of it in accordance with accepted medical practices and the relevant local, state or federal laws and standards.

## 6.4 Magnetic resonance (MR) safety precautions:

In non-clinical studies, it has been shown to be safe to obtain images of Cobalt Chromium L605 stents with ePTFE by Magnetic Resonance (MR) under certain conditions (conditional MR). The MR test conditions used to evaluate these devices were: for interactions of magnetic fields, a static magnetic field with an intensity of 3 Teslas, with a maximum spatial gradient of 2400 gauss/cm; for heating related to MR, a maximum whole-body-averaged specific absorption rate (SAR) of 2.25 W/kg for 15 minutes within a field of 1.5 Teslas and an average specific absorption rate (SAR) of 1.21 W/kg for 15 minutes within a field of 3 Teslas. The increase in temperature (ASTM 2182-11) for a single stent was no more than  $1.98 \pm 0.32^\circ\text{C}$  under the magnetic field of 1.5 Teslas and  $1.16^\circ\text{C} \pm 0.32^\circ\text{C}$  under the magnetic field of 3 Teslas.

The results of the magnetically induced displacement force (ASTM F2052-15) demonstrated a deflection angle of  $7.00^\circ$  and the results of the magnetically induced torque (ASTM F2213-06) did not show any deflection angle at any angular position.

The artefacts in the images (ASTM F21 19-07) were evaluated in a Siemens Skyra 3T MRI scanner. The maximum artefact for the spin-echo sequence was 5.7 mm and for the gradient-echo sequence it was 11.1 mm. The parameters for obtaining images may need to be adjusted for the apparition of artefacts.

## 7. Possible adverse effects/complications

Possible adverse effects and/or complications that may occur before, during or after the procedure include:

- Death
- Myocardial infarction/angina pectoris
- Stroke/embolism/thrombosis
- Emergency bypass surgery
- Endocarditis
- Total occlusion of the artery/thrombosis/embolisation of the stent
- Perforation or dissection of the segment treated
- Air embolism
- Rupture of the balloon
- Renal impairment
- Amputation
- Allergic/anaphylactic reaction to the metals (L605) and/or to the PTFE
- Allergic reaction to the medication, substances or materials used in the intervention (for example anticoagulant or antiplatelet drugs, or to the contrast medium)
- Reocclusion of the area treated: restenosis
- Arterial spasm
- Haemodynamic deterioration
- Oedema
- Fever
- Arrhythmias
- Aneurysm or pseudoaneurysm
- Infection/sepsis
- Hypotension/hypertension
- Local haemorrhage with haematoma at the insertion site
- Arterial thrombus/occlusion near the puncture site
- Arteriovenous fistula
- Abscess
- Tissue reaction/necrosis
- Additional intervention due to, but not limited to:
  - Migration of the endoprosthesis
  - Faulty positioning
  - Partial deployment

## 8. How to use

### 8.1 Materials required

- Heparinised normal saline solution.
- Contrast medium (use contrast media indicated for intravascular use).
- Appropriate medication for the procedure.
- Guidewire of 0.035"; do not use any other size of guidewire.
- Introducer with haemostatic valve of the size indicated on the label. If a long or braided introducer is used, it is advisable to use one French size larger than indicated on the label.
- Three-way valve/extensions.
- Inflation device (manual with pressure gauge).
- Various standard syringes of 10-20 cc with saline solution for cleaning the system.

### 8.2 Preparation of the system

- The stent size must coincide or be slightly larger in diameter than the diameter of the reference vessel and same length than the length of the lesion. Preferably, it should be slightly longer to ensure complete coverage of the lesion.
- Remove the catheter from the protective dispenser. Check that the size is suitable. Remove the protective sheath from the stent.
- Wash the lumen of the guidewire with sterile saline solution injecting it with a syringe from the straight port of the connector until it comes out of the tip. Do not dry with gauze.
- DANGER: If the solution does not exit through the tip of the delivery catheter, the catheter is not in good condition and must not be used.
- **Moisten the stent well with saline solution for 20 seconds.**
- Inspect the stent to check it is located between the radiopaque markers of the balloon. If it is not, it is not use it.
- The internal lumen of the balloon catheter contains air that must be bled before inserting it into the blood stream. If the product is not correctly bled, it can cause the balloon to inflate with air and, consequently, the stent to expand incorrectly. Moreover, in this situation, rupture of the balloon could cause an air embolism in the artery. To bleed the balloon, follow these steps:
  - Fill a syringe or the inflation device with a suitable volume of contrast medium/saline solution. Release the air as indicated by the manufacturer.
  - Connect a three-way valve to an extension line and the extension line to the lateral port of the luer connector of the release system. Connect the syringe to the three-way valve. It is also possible to connect the inflation device to the three-way valve without the extension.
  - Place the catheter on a flat, horizontal surface.
  - With the three-way valve open allowing the flow from the catheter to the syringe/inflation device, induce a negative pressure in the syringe by pulling back on the plunger for a minimum of 30s to let all the air out of the device. The syringe should be in a vertical position.

DANGER: if bubbles do not stop moving from the catheter toward the syringe after 30 seconds at negative pressure, this may be a clear symptom that the balloon catheter has a leak, is broken or the connection between the syringe and three-way valve is not airtight. If you continue to see bubbles after checking the connections, do not use it. Return the device to the manufacturer or distributor for inspection.

- Return the plunger to its initial position, allowing the contrast fluid to penetrate the lumen of the balloon.

DANGER: Do not apply positive pressure as this could cause the balloon to dilate and the stent to expand prematurely.

- Repeat the process again. Once finished, close the three-way valve to stop the passage of air through the balloon.
- Connect the inflation device (if it has not already been connected) filled with 1/2 of a mixture of contrast liquid and saline solution 1:1 (duly bled), to the three-way valve. Do not apply pressure as this could cause the balloon to dilate and the stent to expand prematurely.

### 8.3 Implantation technique

#### Preliminary steps

- Insert the introducer according to the manufacturer's instructions.
- Flush the guidewire with heparinised sterile saline solution. Do not dry with gauze.

- Insert the guidewire through the haemostatic valve of the introducer, in accordance with the manufacturer's instructions. Advance the guidewire inside the vessel using fluoroscopic imaging until it passes the area where the stent will be implanted.

#### Implantation technique

- When the lesion to dilate has a small diameter, is calcified or is so tortuous that it could compromise access to the stent, pre-dilatating it with a balloon catheter must be considered. To pre-dilate the lesion, follow the instructions of the catheter balloon manufacturer. Then implant the stent immediately.
- Insert the proximal end of the guidewire in the distal tip of the release system until the guide passes the straight port of the release system connector.
- Carefully insert the system through the introducer using fluoroscopic imaging. Once the introducer has passed, use fluoroscopic imaging to ensure that the stent has not been displaced from the radiopaque markers of the balloon. The stent has 3 radiopaque markers at each end. These markers must be between the radiopaque markers of the balloon, so that if they are offset, it means that the stent has moved. If this were to happen, the system must be withdrawn in conjunction with the introducer as a single unit in accordance with the instructions in the precautions section of this brochure.
- Carefully advance the stent system on the guidewire, with the aid of fluoroscopy, until you reach the area to be treated. Make sure that the guidewire exits through the tip of the release system.
- Place the stent through the area to be treated; to do this use the radiopaque markers. Before expanding the stent, use high-resolution fluoroscopy to verify that it has not been damaged or moved during the progress. Otherwise, withdraw the system in accordance with the withdrawal procedure indicated in the precautions section of this brochure.
- Slowly inflate the stent release system at low pressure, using the inflation device to apply pressure. Maintain the pressure until you see that the balloon is inflated at the proximal and distal ends of the stent. Continue applying pressure until the nominal pressure indicated in the compliance curve is reached. Do not exceed the rated burst pressure (RBP), which is indicated on the label and compliance curve.
- For optimum expansion, the stent must be in contact with the artery. If it is not, the stent could be displaced when the balloon is removed. The diameter of the balloon must be a little larger than the diameter of the vessel to permit the elastic return of the stent once the balloon has been deflated. Consult the compliance table for a correct expansion as it could cause an arterial dissection if the size is too large.
- Move the plunger of the inflating device back to deflate the balloon. Maintain the negative pressure for 15 to 30 seconds, depending on the size of the balloon. Ensure that the balloon is completely deflated (by fluoroscopy) before moving the catheter.
- With negative pressure in the inflation device and with the guidewire in position, remove the catheter from the lesion and from the patient. Keep the guide through the dilated stenosis.
- Use angiography to confirm the correct placement of the stent. Its radiopaque markers will help indicate that the stent has been positioned correctly.
- Leave the guide running through the dilated stenosis for 30 minutes after angioplasty. If a thrombus is visualised around the stent, administer thrombolytic agents. When the angiogram confirms dilatation, carefully withdraw the guidewire.
- Leave the introducer in place until the haemodynamic profile is back to normal. Suture if necessary, in the usual manner.

#### Additional dilatation of the segments where the product has been implanted

If the size of the stent compared to the vessel is insufficient, the stent postdilatation technique can be applied with its release system or with a balloon with a greater diameter. To do this, insert a balloon catheter in the area of the stent and dilate in accordance with the manufacturer's recommendations.

Stents with diameters of 5 to 8 mm can be dilated after the fact with balloons of up to 10 mm. Stents with diameters of 9 to 10 mm can be dilated after the fact with balloons of up to 12 mm.

#### Withdrawal of the product without expanding

If the stent does not cross the lesion, it may be necessary to withdraw

it unexpanded. To do so, follow the instructions recommended in the precautions of this leaflet and withdraw the system as a single unit.

## Implantation card

An implantation card is included with this system, where the data on the implant, the patient and the hospital must be recorded. To provide the information on the implant, simply place one of the small labels from the label of the bag or packaging. A copy of this card must be submitted to the patient and another copy to the distributor. The patient must carry the card and hand it to any health professional who may attend to them in future.

## 9. Guarantee

The product and all of its components have been designed, manufactured, tested and packaged with maximum precaution. LVD Biotech guarantees the product until its expiry date, provided the package has not been broken, tampered with or damaged.

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CE 0318 © iVascular