



# EFFPAC trial

Multicenter Randomized Controlled Trial to Assess the Effectiveness of Paclitaxel-coated Luminor® Balloon Catheter versus Uncoated Balloon Catheter in the Superficial Femoral and Popliteal Arteries to Prevent Vessel Restenosis or Reocclusion

# Luminor data\_EffPac trial 5-year

---

## Study design

Investigator initiated, prospective, multicenter, randomized controlled trial

## Study objective

To assess efficacy and safety of Luminor-35 paclitaxel-coated balloon angioplasty in SFA/PA lesions up to 60-month (5-year)

## PI


Prof Ulf Teichgräber

## Sponsor

University of Jena, Germany

# Participating Centers

1. Jena	Universitätsklinikum	R. Aschenbach
2. Arnsberg	Klinikum Arnsberg	M. Lichtenberg
3. Bad Krozingen	Herzzentrum	T. Zeller
4. Berlin	Ihre Radiologen	K. Brechtel
5. Hamburg	Angiologikum	S. Sixt
6. Kusel	Westpfalz Klinikum	P. von Flotow
7. Karlsbad	SRH Klinikum	E. Blessing
8. Leipzig	Universitätsklinikum	D. Scheinert
9. München	LMU München	M. Treitl
10. Sonneberg	Medios Kliniken	M. Thieme
11. Heidelberg	Universitätsklinik	B. Vogel



ClinicalTrial.gov Identifier: [NCT02540018](https://clinicaltrials.gov/ct2/show/study/NCT02540018)

# Study Endpoints and Key Eligibility Criteria

## Primary Endpoint

- LLL at 6 months

## Secondary Endpoints

- Binary restenosis
- Primary patency
- Freedom from TLR
- Freedom from TVR
- Rutherford category
- WIQ-score
- ABI
- EQ-5D score

- All-cause mortality
- Target limb amputation


## Inclusion

- Rutherford category 2-4
- De-novo stenotic/restenotic or occluded ( $\geq 70\%$ ) SFA/prox. PA lesions
- Lesion length  $\leq 150$  mm
- 1 lesion/patient
- Successful pre-dilation

## Exclusion

- Previous TV surgery
- Major amputation TL
- Severely calcified lesions (PTA resistant)
- In-stent restenosis

# Patient Flow at 5-year



\*Luminor 35:paclitaxel 3µg/mm<sup>2</sup>

# Baseline Patient Characteristics

	Luminor n = 85	POBA n = 86	p value
Age, years	68.0 ± 7.5	68.1 ± 8.8	p = 0.979
Male, %	60.0	69.8	p = 0.239
Diabetes, %	36.5	40.4	p = 0.681
Hypertension, %	87.1	84.9	p = 0.850
Hyperlipidemia, %	70.7	68.6	p = 0.144
Current smoker, %	40.5	43.0	p = 0.856
Critical limb ischemia, %	3.6	1.2	p = 0.613
ABI	0.73 ± 0.23	0.74 ± 0.23	p = 0.929

# Lesion and Procedure Characteristics

	Luminor (n= 85)	POBA (n= 86)	p value
Lesion length, mm	59.1 ± 43.4	55.8 ± 39.1	p = 0.732
CTO, %	20.2	25.6	p = 0.492
Calcification, %			
Severe	3.6	11.6	p = 0.232
Moderate	42.2	44.2	
Mid / dist. popliteal artery, %	18.8	14.0	p = 0.248
Pre-dilation, %	98.8	98.8	p = 0.993
Dissection, %	37.6	40.7	p = 0.801
Bailout stenting, %	15.3	18.8	p = 0.709
Residual DS, %			
post-angioplasty	15.5 ± 16.7	14.9 ± 16.2	p = 0.807
post-treatment	7.5 ± 9.3	8.3 ± 10.1	p = 0.699


# Primary Endpoint – 6-Month Late lumen loss (LLL)

Study	DCB 6-month LLL	Control 6-month LLL	Difference DCB vs POBA (mm)
<b>THUNDER</b> Tepe et al. 2008 Paccocath coating	0.4±1.2	1.7±1.8	-1.3
<b>AcoArt I</b> Trial Jia et al. 2016 Orchid (Acotec)	0.05±0.73	1.15±0.89	-1.1
<b>EFFPAC 2018</b> Luminor (iVascular)	<b>0.14</b> [CI: -0.38; 0.67]	<b>1.06</b> [CI:0.54; 1.59]	<b>-0.92</b> [CI:-1.364; -0.49] p < 0.001
<b>RANGER</b> Bausback et al. 2017 Ranger DCB	-0.16±0.99	0.76±1.4	-0.92
<b>LEVANT I</b> Scheinert et al. 2014 Lutonix (Bard)	0.46±1.13	1.09±1.07	-0.63
<b>BIOLUX P-I</b> Trial Scheinert et al. 2015 Passeo-18 Lux (Biotronik)	0.51±0.72	1.04±1.0	-0.53
<b>FEMPAC</b> Werk et al. 2008 Paccocath DCB	0.5±1.1	1.0±1.1	-0.5
<b>CONSEQUENT</b> 2017 SeQuent Please (B. Braun)	0.35 [CI: 0.19; 0.79]	0.72 [CI: 0.68; 1.22]	-0.37




# Clinical Improvement: Change of Rutherford Class – 5-year

With Luminor patients clinical improvement is better than with POBA




\*P-value for difference in change from baseline to 60 months between DCB and POBA

# Freedom From TLR \_5-year



Analysing all the existing similar DCB RCT at 5-year follow-up:




Sources:

INPACT SFA: <https://www.ahajournals.org/doi/10.1161/CIRCINTERVENTIONS.118.007702> ; AcoArt I: <https://pubmed.ncbi.nlm.nih.gov/33600928/>; THUNDER: <https://pubmed.ncbi.nlm.nih.gov/25616822/>

# Primary Patency\_5-year

Primary patency after 60M (Kaplan Meier estimate with 95% CI)



**Primary patency:**  
Freedom from restenosis  
(determined by duplex  
ultrasound PSVR < 2.5) and  
freedom from TLR

Same safety outcomes as POBA:

	DCB	POBA	P value
All-cause mortality, %	9* (11.3)	14** (16.3)	p = 0.378
Binary restenosis, %	42 (72.4)	36 (62.1)	p = 0.323
TLR, %	10 (20.4)	18 (35.3)	p = 0.121
Periprocedural complication, %			
Dissection	32 (37.6)	35 (40.7)	p = 0.801
False aneurysm	0	1 (1.2)	p = 1.000
Thromb. embolization	1 (1.2)	0	p = 1.000

\* Reasons: Hyperglycemic coma, stroke, multimorbidity, lung cancer, hear failure, sepsis, cancer, unknown (2x); 5 patients could not be reached

\*\* Reasons: suicide, cardiac arrest (2x), cholangiocellular carcinoma, multiple organ failure, respiratory failure, stroke, pneumonia, heart failure, Covid-19, cardiac insufficiency, unknown (3x)

# EffPac\_ Conclusions

---

- Clinical improvement remains stable over 5 years
- Primary patency with Luminor is superior over POBA throughout 5 years
- No difference in long-term safety between Luminor and POBA was detected
- Luminor angioplasty should be preferred over POBA in femoropopliteal artery disease



One of the few trials with 5-year follow-up

Same safety outcomes than POBA

Significant lesion improvement

**luminor**  
The DCB you need

# Luminor clinical trials

**iVasTriam**  
The iVascular Clinical Trial Program



## LUMINOR registry

Real world registry, 1-year follow-up BTK subgroup. N= 215  
PP: 85.9%; fTLR:89.6%

## EffPac

Randomized controlled trial (Luminor vs PTA). N=171  
5-year follow-up. fTLR=82.1%, PP=61.4%

## TINTIN

Prospective trial in complex SFA lesions. N=100  
3-year follow-up. fTLR= 82.2%

## MERLION

Prospective BTK trial. N= 50  
1-year follow-up. fTLR= 81.6%

## BIBLIOS

Prospective BTK trial, N=150  
Preliminary 1-year follow-up. fTLR= 88%

## LUMIFOLLOW

France Registry. N=500

## SOL-JAPAN

Prospective, single arm trial in Japan. N=120



---

# Thank you