

Cera™

Vascular Plug System



Fast and Complete Occlusion

The Cera™ Vascular Plug is indicated for arterial and venous embolizations in the peripheral vasculature.



CE 0344

FEATURES & BENEFITS

Great Biocompatibility

- Titanium Nitride (TiN) coating prevents thrombosis and accelerates endothelialization

Fast and Complete Occlusion

- PTFE membrane reduces total occlusion time and eliminates residual shunt


Easy and Precise Deployment

- Stainless steel microscrew connection

Available in Various Sizes Ranging

- Allows for use in various vessel lesions

ORDERING INFORMATION



Code	Plug Diameter (D) (mm)	Unconstrained Plug Length (L) (mm)	Minimum ID per Cera™ VP size (Fr)
LT-PLUG-04	4	7	4F(1.40mm/0.0551")
LT-PLUG-06	6	7	4F(1.40mm/0.0551")
LT-PLUG-08	8	7	5F(1.75mm/0.0689")
LT-PLUG-10	10	8	5F(1.75mm/0.0689")
LT-PLUG-12	12	8	6F(2.10mm/0.0827")
LT-PLUG-14	14	10	7F(2.40mm/0.0945")
LT-PLUG-16	16	10	7F(2.40mm/0.0945")
LT-PLUG-18	18	11	8F(2.80mm/0.110")
LT-PLUG-20	20	12	8F(2.80mm/0.110")
LT-PLUG-22	22	14	9F(3.10mm/0.122")
LT-PLUG-24	24	14	9F(3.10mm/0.122")

Lifetech Scientific (Shenzhen) Co., Ltd.

Add: Cybio Electronic Building, Langshan 2nd Street, North Area of High-tech Park, Nanshan District, Shenzhen 518057, P. R. China

Tel: +86-755-86026250 Fax: +86-755-86026251

Email: marketing@lifetechmed.com

Web: www.lifetechmed.com

EU authorized representative:
LifeTech Scientific (Europe) Coöperatief U.A.
Kruisdonk 64, 6222 PH Maastricht, the Netherlands