

**iVascular**<sup>®</sup>  
therapies for living

# navitian

Coronary microcatheter



The navigation you need



[www.ivascular.global](http://www.ivascular.global)

Navitian, the needed microcatheter to arrive to the most challenging lesions.

**Its exceptional design brings benefits as:**

**Outstanding trackability**  
through the small and  
tortuous arteries due to  
hydrophilic coating  
HYDRAX PLUS


**Optimal flexibility,  
pushability and resistance  
to kinking** provided by a  
full length variable braided  
pattern

**Complete  
guidewire control**  
ensured by its  
conical transition

**Superior crossability and  
penetration capacity** due to  
its optimized tip profile and  
shape with double c-cut  
radiopaque markers

## navitian features

- > **Over the wire catheter (OTW)**
- > **Guidewire compatibility** 0.014"
- > **Usable catheter length** 135 and 150 cm
- > **Radiopaque markers**
- > **Internal PTFE layer**
- > **Hydrophilic coating:** 70 and 90 cm



Product with CE mark, certified by Notified Body 0318


### References

Usable length 135 cm
Usable length 150 cm

MCC C14 135 001
MCC C14 150 001

The availability of each reference for the sale is linked to the authorization of commercialization in the country of destination

Distributed by:



Manufactured by:

**Life Vascular Devices Biotech S.L**  
[www.ivascular.global](http://www.ivascular.global)  
[info@ivascular.global](mailto:info@ivascular.global)